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1. Abrahamson M et al: Biochem J 1990;268:287-94.
2. Laterza OF et al: Clin Chem 2002;48:63-99.
3. Grubb AO. Adv Clin Chem 2000;35:63-99.
4. Filler G et al: Clin Biochem 2005;38:1-8.
5. Flodin M et al: Scand J Clin Lab Invest 2007;67:560-567.
6. Shlipak M.G, et al: Clinical Chemistry 57: 737-745, 2011.
7. CLSI. Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline – Second Edition. CLSI document EP17-A2. Wayne, PA: Clinical and Laboratory Standards Institute; 2012.
8. CLSI. Evaluation of the Linearity of Quantitative Measurement Procedures. 2nd ed. CLSI guideline EP06. Clinical Laboratory Standards Institute; 2020.
9. CLSI. Establishing and Verifying and Extended Measuring Interval Through Specimen Dilution and Spiking; 1st Edition. CLSI guideline EP34. Wayne, PA: Clinical and Laboratory Institute; 2018.
10. CLSI. Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline – Third Edition. CLSI document EP05-A3. Wayne, PA: Clinical Laboratory Standards Institute; 2014.
11. Westgard JO. Basic Method Validation, 3rd Edition. 2008; ISBN13: 9781886958258.
12. Sonntag O, Scholer A. Ann Clin Biochem 2001;38:376-85.
13. Larsson A et al: J Immunol Methods. 1988 Apr 6;108(1-2):205-8.
14. CLSI. Interference Testing in Clinical Chemistry. 3rd ed. CLSI guideline EP07. Wayne, PA: Clinical Laboratory Standards Institute; 2018.
15. CLSI. Measurement Procedure Comparison and Bias Estimation Using Patient Samples. 3rd ed. CLSI guideline EP09c. Wayne, PA: Clinical Laboratory Standards Institute; 2018.
16. Qiu X et al: Oncotarget. 2017;8(42):72985-72999.
17. CLSI. Defining, Establishing, and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline – Third Edition. CLSI document C28-A3c. Wayne, PA: Clinical Laboratory Standards Institute; 2008.

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Instrument Settings for Cystatin C Immunoassay

Cystatin C AU5800 application settings

System Reagent: B08179 Reagent ID: 228

Parameters		Specific Test Parameters			
General	LIH	ISE	Calculated Test	Range	
Test Name: <input type="text" value="CysC"/> < > Type: <input type="text" value="Serum"/> Operation <input type="text" value="Yes"/>					
Sample Volume	<input type="text" value="2"/> μL	Dilution	<input type="text" value="0"/> μL	OD Limit	
Pre-Dilution Rate	<input type="text" value="1"/>	Diluent Bottle	<input type="text" value="Outside"/>	Min.OD	<input type="text"/>
Rgt. Volume	R1(R1-1) <input type="text" value="150"/> μL	Dilution	<input type="text" value="0"/> μL	Reagent OD Limit	
	R1-2 <input type="text"/>	Dilution	<input type="text"/>	First Low	<input type="text" value="-2.0"/> High <input type="text" value="2.0"/>
	R2(R2-1) <input type="text" value="30"/> μL	Dilution	<input type="text" value="10"/> μL	Last Low	<input type="text"/>
Common Rgt. Type	<input type="text"/>	Name	<input type="text"/>	Dynamic Range Low	<input type="text" value="0.49"/> High <input type="text" value="7.07"/>
Wavelength	Pri <input type="text" value="540"/> nm	Sec.	<input type="text"/>	Correlation Factor A	<input type="text" value="1.00"/> B <input type="text" value="0.00"/>
Method	<input type="text" value="End Point"/>			Factor for Maker A	<input type="text"/>
Reaction Slope	<input type="text" value="+"/>			Onboard Stability Period	<input type="text" value="60**"/> Day <input type="text"/> Hour
Measuring Point1 First	<input type="text" value="13"/>	Last	<input type="text" value="27"/>	LIH Influence Check	<input type="text"/>
Measuring Point2 First	<input type="text"/>	Last	<input type="text"/>	Lipemia	<input type="text"/>
Linearity Limit	<input type="text"/>			Icterus	<input type="text"/>
Lag Time Check	<input type="text"/>			Hemolysis	<input type="text"/>

Parameters		Specific Test Parameters					
General	LIH	ISE	Calculated Test	Range			
Test Name: <input type="text" value="CysC"/> < > Type: <input type="text" value="Serum"/>							
Value/Flag: <input type="text" value="#"/>							
Level Low <input type="text" value="#"/> High <input type="text" value="#"/>							
Specific Ranges:							
	Sex	Year	Month	Year	Month	Low	High
<input type="checkbox"/> 1.	<input type="text" value="#"/>						
<input type="checkbox"/> 2.	<input type="text" value="#"/>						
<input type="checkbox"/> 3.	<input type="text" value="#"/>						
<input type="checkbox"/> 4.	<input type="text" value="#"/>						
<input type="checkbox"/> 5.	<input type="text" value="#"/>						
<input type="checkbox"/> 6.	<input type="text" value="#"/>						
<input type="checkbox"/> 7.	Standard demographics						
<input type="checkbox"/> 8.	Not within expected values						
Panic Value		Low	<input type="text"/>	High	<input type="text"/>	Unit	<input type="text" value="mg/L"/>
						Decimal Places	<input type="text" value="#"/>

Parameters		Calibration Parameters			
Calibrators	Calibration Specific	STAT Table Calibration			
General	ISE				
Test Name: <input type="text" value="CysC"/> < > Type: <input type="text" value="Serum"/> Cuvette: <input type="text"/>					
<input type="checkbox"/> Use Serum Cal.					
Calibration Type: <input type="text" value="GAB"/>		Formula: <input type="text" value="Spline"/>		Counts: <input type="text" value="#"/>	
<Calibrator Parameters>					
	Calibrator	OD	Conc	Factor Range	
				Low	High
Point 1:	<input type="text" value="1"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
Point 2:	<input type="text" value="2"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
Point 3:	<input type="text" value="3"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
Point 4:	<input type="text" value="4"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
Point 5:	<input type="text" value="5"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
Point 6:	<input type="text" value="6"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
Point 7:	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
Point 8:	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
Point 9:	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
Point 10:	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
<Point Cal. For Master Curve>		No. of Correction Points	<input type="text"/>	<input type="checkbox"/> Use Master Curve	<input type="checkbox"/>
		OD Range	<input type="text"/>	<input type="checkbox"/> Lot Calibration	<input type="checkbox"/>
	Calibrator	OD	Conc	Low	High
Point-1	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
Point-2	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
MB Type Factor: <input type="text"/>		1-Point Calibration Point <input type="text"/>		<input type="checkbox"/> with Conc-0	
				Slope Check	<input type="text" value="+"/>
				Allowance Range Check	
				<input type="checkbox"/> Reagent Blank	<input type="text"/>
				<input type="checkbox"/> Calibration	<input type="text"/>
				Advanced Calibration Operation	<input type="text" value="#"/>
				Interval (RB/ACAL)	<input type="text" value="#"/>
				Stability	
				Reagent Blank	<input type="text" value="28"/> Day <input type="text" value="0"/> Hour
				Calibration	<input type="text" value="28"/> Day <input type="text" value="0"/> Hour

User defined
 * Lot specific, see analytical value sheet available on www.gentian.com
 ** Based on results from instrument AU400 (Beckman Coulter®)

Cystatin C AU680 application settings
System Reagent: B08179 Reagent ID: 228

Parameters		Specific Test Parameters						
General	LIH	ISE	Calculated Test	Range				
Test Name:	CysC ▾	<	>	Type:	Serum ▾	Operation	Yes ▾	
Sample Volume	2 μL	Dilution	0 μL	OD Limit				
Pre-Dilution Rate	1 ▾			Min.OD		Max.OD		
Rgt. Volume R1(R1-1)	150 μL	Dilution	0 μL	Reagent OD Limit				
				First	Low	-2.0	High	2.0
				Last	Low		High	
R2(R2-1)	30 μL	Dilution	10 μL	Dynamic Range Low	0.44	High	7.30	
Common Rgt. Type		Name		Correlation Factor A	1.00	B	0.00	
Wavelength	Pri 540 ▾nm	Sec.		Factor for Maker A		B		
Method	End Point ▾							
Reaction Slope	+ ▾			Onboard Stability Period	60**	Day		
Measuring Point1 First	13	Last	27	LIH Influence Check				
Measuring Point2 First		Last		Lipemia				
Linearity Limit				Icterus				
Lag Time Check				Hemolysis				

Parameters		Specific Test Parameters						
General	LIH	ISE	Calculated Test	Range				
Test Name:	CysC ▾	<	>	Type:	Serum ▾			
Value/Flag:	# ▾	Low	#	High	#			
Level	#	From	#	To	#			
Specific Ranges:		Sex	Year	Month	Year	Month	Low	High
<input type="checkbox"/> 1.	# ▾	# ▾	# ▾	# ▾	# ▾	# ▾	# ▾	# ▾
<input type="checkbox"/> 2.	# ▾	# ▾	# ▾	# ▾	# ▾	# ▾	# ▾	# ▾
<input type="checkbox"/> 3.	# ▾	# ▾	# ▾	# ▾	# ▾	# ▾	# ▾	# ▾
<input type="checkbox"/> 4.	# ▾	# ▾	# ▾	# ▾	# ▾	# ▾	# ▾	# ▾
<input type="checkbox"/> 5.	# ▾	# ▾	# ▾	# ▾	# ▾	# ▾	# ▾	# ▾
<input type="checkbox"/> 6.	# ▾	# ▾	# ▾	# ▾	# ▾	# ▾	# ▾	# ▾
<input type="checkbox"/> 7.	No demographics							
<input type="checkbox"/> 8.	Not within expected values							
Unit	mg/L	Decimal Places	#	Panic Value				
				Low	#	High	#	

Parameters		Calibration Parameters					
General	Calibration Specific	STAT Table Calibration					
Test Name:	CysC ▾	<	>	Type	Serum ▾	<input type="checkbox"/> Use Serum Cal.	
Calibration Type:	6AB ▾	Formula:	Spline ▾	Counts:	# ▾		
<Calibrator Parameters>							
Calibrator	OD	Conc	Low	High	Slope Check	+ ▾	
Point 1:	1 ▾				Allowance Range Check		
Point 2:	2 ▾				<input type="checkbox"/> Reagent Blank		
Point 3:	3 ▾				<input type="checkbox"/> Calibration		
Point 4:	4 ▾				Advanced Calibration		
Point 5:	5 ▾				Operation	# ▾	
Point 6:	6 ▾				Interval (RB/ACAL)	# ▾	
Point 7:	▾				<input type="checkbox"/> Use Master Curve	<input type="checkbox"/> Lot Calibration	
Point 8:	▾				Master Curve>		
Point 9:	▾				No. of Correction Points	▾	
Point 10:	▾				OD Range		
Point-1	▾	OD	Conc	Low	High	Stability	
Point-2	▾					Reagent Blank	28 Day ▾ Hour
						Calibration	28 Day ▾ Hour
MB Type Factor:		1-Point Calibration Point	▾	<input type="checkbox"/> with Conc-0			

User defined

* Lot specific, see analytical value sheet available on www.gentian.com

** Based on results from instrument AU400 (Beckman Coulter®)

Cystatin C AU480 application settings
System Reagent: B08179 Reagent ID: 228

Parameters		Specific Test Parameters			
General	LIH	ISE	Calculated Test	Range	
Test Name: <input type="text" value="CysC"/> < > Type: <input type="text" value="Serum"/> Operation: <input type="text" value="Yes"/>					
Sample Volume	<input type="text" value="2"/> μ L	Dilution	<input type="text" value="0"/> μ L	OD Limit	
Pre-Dilution Rate	<input type="text" value="1"/>			Min.OD	<input type="text"/>
Rgt. Volume	R1(R1-1) <input type="text" value="150"/> μ L	Dilution	<input type="text" value="0"/> μ L	Reagent OD Limit	
				First Low	<input type="text" value="-2.0"/> High <input type="text" value="2.0"/>
				Last Low	<input type="text"/>
		R2(R2-1) <input type="text" value="30"/> μ L	Dilution <input type="text" value="10"/> μ L	Dynamic Range Low	<input type="text" value="0.43"/> High <input type="text" value="7.32"/>
Common Rgt. Type		Name	<input type="text"/>	Correlation Factor A	<input type="text" value="1.00"/> B <input type="text" value="0.00"/>
Wavelength	Pri <input type="text" value="540"/> nm	Sec.	<input type="text"/>	Factor for Maker A	<input type="text"/>
Method	<input type="text" value="End Point"/>				
Reaction Slope	<input type="text" value="+"/>			Onboard Stability Period	<input type="text" value="60**"/> Day <input type="text"/> Hour
Measuring Point1 First	<input type="text" value="13"/>	Last	<input type="text" value="27"/>	LIH Influence Check	<input type="text"/>
Measuring Point2 First	<input type="text"/>	Last	<input type="text"/>	Lipemia	<input type="text"/>
Linearity Limit	<input type="text"/>			Icterus	<input type="text"/>
Lag Time Check	<input type="text"/>			Hemolysis	<input type="text"/>

Parameters		Specific Test Parameters					
General	LIH	ISE	Calculated Test	Range			
Test Name: <input type="text" value="CysC"/> < > Type: <input type="text" value="Serum"/>							
Value/Flag:	<input type="text" value="#"/>	Low	<input type="text"/>	High	<input type="text"/>		
Level							
Specific Ranges:	From	To		Panic Value			
				Low	High		
<input type="checkbox"/> 1.	Sex <input type="text" value="#"/>	Year <input type="text" value="#"/>	Month <input type="text" value="#"/>	Year <input type="text" value="#"/>	Month <input type="text" value="#"/>	Low <input type="text" value="#"/>	High <input type="text" value="#"/>
<input type="checkbox"/> 2.	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
<input type="checkbox"/> 3.	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
<input type="checkbox"/> 4.	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
<input type="checkbox"/> 5.	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
<input type="checkbox"/> 6.	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
<input type="checkbox"/> 7.	No demographics					<input type="text" value="#"/>	<input type="text" value="#"/>
<input type="checkbox"/> 8.	Not within expected values					<input type="text" value="#"/>	<input type="text" value="#"/>
Unit	<input type="text" value="mg/L"/>	Decimal Places	<input type="text" value="#"/>				

Parameters		Calibration Parameters			
Calibrators	Calibration Specific	STAT Table Calibration			
General	ISE				
Test Name: <input type="text" value="CysC"/> < > Type: <input type="text" value="Serum"/> <input type="checkbox"/> Use Serum Cal.					
Calibration Type: <input type="text" value="6AB"/>		Formula: <input type="text" value="Spline"/>		Counts: <input type="text" value="#"/>	
<Calibrator Parameters>					
Calibrator	OD	Conc	Low	High	Slope Check <input type="text" value="+"/>
Point 1:	<input type="text" value="1"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	Allowance Range Check
Point 2:	<input type="text" value="2"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="checkbox"/> Reagent Blank <input type="text"/>
Point 3:	<input type="text" value="3"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="checkbox"/> Calibration <input type="text"/>
Point 4:	<input type="text" value="4"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	Advanced Calibration
Point 5:	<input type="text" value="5"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	Operation <input type="text" value="#"/>
Point 6:	<input type="text" value="6"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	Interval (RB/ACAL) <input type="text" value="#"/>
Point 7:	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	
Point 8:	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	
Point 9:	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	
Point 10:	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	
<Point Cal. For	No. of Correction Points <input type="text"/>	<input type="checkbox"/> Use Master Curve	<input type="checkbox"/>	<input type="checkbox"/> Lot Calibration	
<Master Curve>					
Calibrator	OD	Conc	Low	High	Stability
Point-1	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	Reagent Blank <input type="text" value="28"/> Day <input type="text"/>
Point-2	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	Calibration <input type="text" value="28"/> Day <input type="text"/>
MB Type Factor:	<input type="text"/>	1-Point Calibration Point <input type="checkbox"/>	<input type="text"/>	<input type="checkbox"/> with Conc-0	

User defined
 * Lot specific, see analytical value sheet available on www.gentian.com
 ** Based on results from instrument AU400 (Beckman Coulter®)

Cystatin C Dx C 500 AU application settings System Reagent: B08179 Reagent ID: 228

TEST CONFIGURATION & CHEMISTRY DETAILS												
Assay Name	Test	Rev	Discipline		Chemistry							
Test ID	CYS		Calculated Result		<input type="checkbox"/>							
LIS Code	CYS											
UNITS AND RANGE SETTINGS												
Use Settings from	Serum	Units	mg/L	Decimal Places	x.xx	Plasma						
Test Kind	General	Revision	01	<input checked="" type="checkbox"/> Multi Reagent Switch								
Reagent Name	CYS	Reagent ID	228	<input type="checkbox"/> FSE Test								
ABB Name	CYS1G	Parameter Long Name										
Cystatin C B08179 CYS1G CYS C Serum												
Region	<input checked="" type="checkbox"/> US <input checked="" type="checkbox"/> OUS <input checked="" type="checkbox"/> AP <input type="checkbox"/> JP <input checked="" type="checkbox"/> EU <input type="checkbox"/> Other											
GENERAL PARAMETERS												
SAMPLE VOLUME		Sample Volume	2.0	µL	Dilution	0	µL	REACTION OD LIMIT				
REAGENT VOLUME		Predilution Rate	1				Low		High			
WAVELENGTH		R1-1	150	µL	Dilution	0	µL	REACTION BLANK OD LIMIT				
METHOD		R2-1	30	µL	Dilution	10	µL	First Low	-2.0000	High	2.0000	
REACTION SLOPE								Last Low	-2.0000	High	2.0000	
MEASURING POINT								ANALYTICAL MEASURING RANGE				
Linearity Limit								Low	0.38	High	7.84	
Lag Time Check								MANUFACTURER FACTOR				
								A	1	B	0	
								REAGENT ONBOARD STABILITY				
								60** Days			0** Hours	
								LIH INFLUENCE CHECK				
								<input type="checkbox"/> Perform LIH check				
								Lipemia			+ ▼	
								Icterus			+ ▼	
								Hemolysis			+ ▼	
CALIBRATION PARAMETERS												
Base Unit	Decimal Place	Unit 1	Factor 1	Unit 2	Factor 2	Unit 3	Factor 3	Unit 4	Factor 4			
mg/L	2	None	0	None	0	None	0	None	0			
CALIBRATOR SPECIFIC												
Calibration Type		6AB		Counts	2		<input type="checkbox"/> Use highest calibrator for Upper AMR					
Formula		Spline		MB Factor			CALIBRATION OD AND CONCENTRATION PARAMETERS					
Calibrator Name		Add CYS		Positive Cutoff			Point 1	Calibrator Name	Conc	OD Range Low	OD Range High	
<input checked="" type="checkbox"/> SLOPE CHECK		Number of Levels 6						Point 2	CYS CAL-2	*	-2.0000	2.0000
Slope Check		+						Point 3	CYS CAL-3	*	-2.0000	2.0000
STABILITY AND INTERVAL								Point 4	CYS CAL-4	*	-2.0000	2.0000
Reagent Blank Stability		28	Days	0	Hours	Interval	Lot	Point 5	CYS CAL-5	*	-2.0000	2.0000
Calibration Stability		28	Days	0	Hours	Interval	Lot	Point 6	CYS CAL-6	*	-2.0000	2.0000
								Point 7				
								OD DELTA CHECK				
								<input type="checkbox"/> Reagent Blank				0.0000
								<input type="checkbox"/> Calibration				0.0000
PROZONE CHECK PARAMETERS												
<input type="checkbox"/> Logio Check 1								<input type="checkbox"/> Logio Check 2				
Check Points	Point 1	#	Decision Values	Value 1	#	Point 1	#	Check Points	Point 1	#	Decision Values	
	Point 2	#	Value 2	#	Interval	#	Value 1	#	Interval	#	Value 1	
	Point 3	#	Value 3	#			Value 2	#			Value 2	
Limit Points		Limit 1	#	Limit Points		Limit 1	#	Limit Points		Limit 1	#	
		Limit 2	#			Limit 2	#			Limit 2	#	
Check Pattern		Pattern										

User defined
* Lot specific, see analytical value sheet available on www.gentian.com
** Based on results from instrument AU400 (Beckman Coulter®)

Cystatin C DxC 700 AU application settings
System Reagent: B08179 Reagent ID: 228

General	LIH	ISE	Calculated Test	Range
Test Name: <input type="text" value="CYS1G"/> <input type="button" value="Test No"/> Type: <input type="text" value="Serum"/> Operation: <input type="text" value="Yes"/>				
Sample Volume <input type="text" value="2.0"/> μ L Pre-Dilution Rate <input type="text" value="1"/>	Dilution <input type="text" value="0"/> μ L	OD Limit	Min OD <input type="text"/> Max OD <input type="text"/>	
Reagent Volume R1 (R1-1) <input type="text" value="150"/> μ L R1-2 <input type="text"/> μ L R2 (R2-1) <input type="text" value="30"/> μ L	Dilution <input type="text" value="0"/> μ L Dilution <input type="text"/> μ L Dilution <input type="text" value="10"/> μ L	Reagent OD Limit 1* Low <input type="text" value="-2.0000"/> High <input type="text" value="2.0000"/> Last Low <input type="text" value="-2.0000"/> High <input type="text" value="2.0000"/> Analytical Measuring Range Low <input type="text" value="0.40"/> High <input type="text" value="8.07"/>		
Common Reagent Type <input type="text" value="None"/> Name <input type="text" value="None"/>		Correlation Factor A <input type="text" value="1"/> B <input type="text" value="0"/>		
Wavelength Pri <input type="text" value="540"/> nm Method <input type="text" value="END"/>	Sec <input type="text" value="None"/> nm	Manufacturer Factor A <input type="text" value="1"/> B <input type="text" value="0"/>		
Reaction Slope <input type="text" value="+"/>		Onboard Stability Period <input type="text" value="60"/> Day <input type="text" value="0"/> Hour		
Measuring Point-1 1st <input type="text" value="13"/> Last <input type="text" value="27"/>		LIH Influence Check <input type="text" value="No"/>		
Measuring Point-2 1st <input type="text"/> Last <input type="text"/>		Lipemia <input type="text"/>		
Linearity Limit <input type="text"/> %		Icterus <input type="text"/>		
Lag Time Check <input type="text"/>		Hemolysis <input type="text"/>		

General	LIH	ISE	Calculated Test	Range	
Test Name: <input type="text" value="CYS1G"/> <input type="button" value="Test No"/> Type: <input type="text" value="Serum"/>					
Value/Flag <input type="text" value="Value"/>		Level Low <input type="text" value="99999.99"/> High <input type="text" value="99999.99"/>			
Specific Ranges					
<input type="checkbox"/> 1:	Sex <input type="text" value="#"/>	From Year <input type="text" value="#"/> Month <input type="text" value="#"/>	To Year <input type="text" value="#"/> Month <input type="text" value="#"/>	Other Type <input type="text" value="None"/>	Low <input type="text" value="#"/> High <input type="text" value="#"/>
<input type="checkbox"/> 2:	<input type="text" value="#"/>	<input type="text" value="#"/> <input type="text" value="#"/>	<input type="text" value="#"/> <input type="text" value="#"/>	<input type="text" value="None"/>	<input type="text" value="#"/> <input type="text" value="#"/>
<input type="checkbox"/> 3:	<input type="text" value="#"/>	<input type="text" value="#"/> <input type="text" value="#"/>	<input type="text" value="#"/> <input type="text" value="#"/>	<input type="text" value="None"/>	<input type="text" value="#"/> <input type="text" value="#"/>
<input type="checkbox"/> 4:	<input type="text" value="#"/>	<input type="text" value="#"/> <input type="text" value="#"/>	<input type="text" value="#"/> <input type="text" value="#"/>	<input type="text" value="None"/>	<input type="text" value="#"/> <input type="text" value="#"/>
<input type="checkbox"/> 5:	<input type="text" value="#"/>	<input type="text" value="#"/> <input type="text" value="#"/>	<input type="text" value="#"/> <input type="text" value="#"/>	<input type="text" value="None"/>	<input type="text" value="#"/> <input type="text" value="#"/>
<input type="checkbox"/> 6:	<input type="text" value="#"/>	<input type="text" value="#"/> <input type="text" value="#"/>	<input type="text" value="#"/> <input type="text" value="#"/>	<input type="text" value="None"/>	<input type="text" value="#"/> <input type="text" value="#"/>
7:	Standard demographics				<input type="text" value="#"/> <input type="text" value="#"/>
8:	Not within expected values				<input type="text" value="#"/> <input type="text" value="#"/>
Critical Limits Low <input type="text" value="#"/> High <input type="text" value="#"/>		Unit <input type="text" value="mg/L"/>	<input type="button" value="Select"/>	Decimal Places <input type="text" value="2"/>	

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