Gentian Cystatin C Immunoassay na sistemima AU (AU5800, AU680, AU480 , DxC 500 AU/DxC 500i, DxC 700 AU) kompanije Beckman Coulter®

REF B08179

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2025-04-03

Za druge jezike posetite:

www.gentian.com/products/ifu/cystatin-c/beckmancoulter

Predstavnici

Odgovorna osoba u UK
Emergo Consulting (UK) Limited
c/o Cr360 – UL International
Compass House, Vision Park Histon
Cambridge CB24 9BZ
Ujedinjeno Kraljevstvo

MedEnvoy Switzerland
Gothardstrasse 28
6302 Zug
Švajcarska

Literatura

1. Abrahamson M et al: Biochem J 1990;268:287-94.
2. Laterza OF et al: Clin Chem 2002;48:63-99.
3. Grubb AO. Adv Clin Chem 2000;35:63-99.
4. Filler G et al: Clin Biochem 2005 ;38 :1-8.
5. Flodin M et al: Scand J Clin Lab Invest 2007;67:560-567.
6. Shlipak M.G, et al: Clinical Chemistry 57: 737-745, 2011.
7. CLSI. *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline – Second Edition*. CLSI document EP17-A2. Wayne, PA: Clinical and Laboratory Standards Institute ;2012.
8. CLSI. *Evaluation of the Linearity of Quantitative Measurement Procedures*. 2nd ed. CLSI guideline EP06. Clinical Laboratory Standards Institut;2020.
9. CLSI. *Establishing and Verifying and Extended Measuring Interval Through Specimen Dilution and Spiking*; 1st Edition. CLSI guideline EP34. Wayne, PA: Clinical and Laboratory Institute; 2018.
10. CLSI. *Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline – Third Edition*. CLSI document EP05-A3. Wayne, PA: Clinical Laboratory Standards Institute; 2014.
11. Westgard JO. Basic Method Validation, 3rd Edition. 2008; ISBN13: 9781886958258.
12. Sonntag O, Scholer A. Ann Clin Biochem 2001;38:376-85.
13. Larsson A et al: J Immunol Methods. 1988 Apr 6;108(1-2):205-8.
14. CLSI. *Interference Testing in Clinical Chemistry*. 3rd ed. CLSI guideline EP07. Wayne, PA: Clinical Laboratory Standards Institute; 2018.
15. CLSI. *Measurement Procedure Comparison and Bias Estimation Using Patient Samples*. 3rd ed. CLSI guideline EP09c. Wayne, PA: Clinical Laboratory Standards Institute; 2018.
16. Qiu X et al.: Oncotarget. 2017;8(42):72985-72999.
17. CLSI. *Defining, Establishing, and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline – Third Edition*. CLSI document C28-A3c. Wayne, PA: Clinical Laboratory Standards Institute; 2008.
18. Nitsch D, et al. Am J Kidney Dis. Jun 2011;57(6):863-72.

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Sistemski reagens: B08179 ID reagensa: 228

General	LIH	ISE	Calculated Test	Range
Test Name: <input type="text" value="CYS1G"/> <input type="text" value="Test No"/> Type: <input type="text" value="Serum"/> Operation: <input type="text" value="Yes"/>				
Sample Volume	<input type="text" value="2.0"/> μL	Dilution	<input type="text" value="0"/> μL	OD Limit
Pre-Dilution Rate	<input type="text" value="1"/>			Min. OD <input type="text"/> Max OD <input type="text"/>
Reagent Volume	R1 (R1-1) <input type="text" value="150"/> μL	Dilution	<input type="text" value="0"/> μL	Reagent OD Limit ^{1st} Low <input type="text" value="-2.0000"/> High <input type="text" value="2.0000"/>
	R1-2 <input type="text"/> μL	Dilution	<input type="text"/> μL	Last Low <input type="text" value="-2.0000"/> High <input type="text" value="2.0000"/>
	R2 (R2-1) <input type="text" value="30"/> μL	Dilution	<input type="text" value="10"/> μL	Analytical Measuring Range Low <input type="text" value="0.40"/> High <input type="text" value="8.07"/>
Common Reagent	Type <input type="text" value="None"/>	Name	<input type="text" value="None"/>	Correlation Factor A <input type="text" value="1"/> B <input type="text" value="0"/>
Wavelength	Pri <input type="text" value="540"/> nm	Sec	<input type="text" value="None"/> nm	Manufacturer Factor A <input type="text" value="1"/> B <input type="text" value="0"/>
Method	<input type="text" value="END"/>			Onboard Stability Period <input type="text" value="60"/> Day <input type="text" value="0"/> Hour
Reaction Slope	<input type="text" value="+"/> μL			LIH Influence Check <input type="text" value="No"/>
Measuring Point-1	1st <input type="text" value="13"/>	Last	<input type="text" value="27"/>	Lipemia <input type="text"/>
Measuring Point-2	1st <input type="text"/>	Last	<input type="text"/>	Icterus <input type="text"/>
Linearity Limit	<input type="text"/> %			Hemolysis <input type="text"/>
Lag Time Check	<input type="text"/>			

General	LIH	ISE	Calculated Test	Range
Test Name: <input type="text" value="CYS1G"/> <input type="text" value="Test No"/> Type: <input type="text" value="Serum"/>				
Value/Flag	<input type="text" value="Value"/>	Level	Low <input type="text" value="-99999.99"/> High <input type="text" value="99999.99"/>	
Specific Ranges				
	Sex <input type="text" value="#"/>	Year From <input type="text" value="#"/>	Month <input type="text" value="#"/>	Year To <input type="text" value="#"/>
	Month <input type="text" value="#"/>	Year <input type="text" value="#"/>	Month <input type="text" value="#"/>	Other Type <input type="text" value="None"/>
<input type="checkbox"/> 1:	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
<input type="checkbox"/> 2:	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
<input type="checkbox"/> 3:	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
<input type="checkbox"/> 4:	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
<input type="checkbox"/> 5:	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
<input type="checkbox"/> 6:	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
7:	Standard demographics			<input type="text" value="#"/>
8:	Not within expected values			<input type="text" value="#"/>
Critical Limits	Low <input type="text" value="#"/>	High <input type="text" value="#"/>	Unit <input type="text" value="mg/L"/>	Select <input type="text" value="2"/> Decimal Places

Calibrators	General	ISE		
Test Name: <input type="text" value="CYS1G"/> <input type="checkbox"/> Use Serum Cal. Type: <input type="text" value="Serum"/>				
Calibration Type:	<input type="text" value="6AB"/>	Formula: <input type="text" value="Spline"/> Counts: <input type="text" value="2"/>		
<Calibrator Parameters>				
	Calibrator	OD		
	Conc	Range		
		Low		
		High		
Point-1	CYSC Calibrator Level 1	*	-2.0000	2.0000
Point-2	CYSC Calibrator Level 2	*	-2.0000	2.0000
Point-3	CYSC Calibrator Level 3	*	-2.0000	2.0000
Point-4	CYSC Calibrator Level 4	*	-2.0000	2.0000
Point-5	CYSC Calibrator Level 5	*	-2.0000	2.0000
Point-6	CYSC Calibrator Level 6	*	-2.0000	2.0000
Point-7				
MB Type Factor	<input type="text"/>	1-Point Calibration Point <input type="text" value="None"/>	<input type="checkbox"/> with Conc-0	Slope Check <input type="text" value="+"/> Allowable Range Check <input type="checkbox"/>
				Reagent Blank <input type="text"/>
				Calibration <input type="text"/>
				Advanced Calibration <input type="text"/>
				Operation <input type="text" value="No"/>
				Interval (RB) <input type="text"/>
				Interval (ACAL) <input type="text"/>
				Stability
				Reagent Blank <input type="text" value="28"/> Day <input type="text" value="0"/> Hour
				Calibration <input type="text" value="28"/> Day <input type="text" value="0"/> Hour

Definiše korisnik

* Specifično za seriju, pogledajte tabelu sa analitičkom vrednošću koja je dostupna na www.gentian.com

** Na osnovu rezultata za instrument AU400 (Beckman Coulter®)