

Cystatin C



Gentian Cystatin C Immunoassay on Beckman Coulter® AU Systems (AU5800, AU680, AU480, DxC 500 AU, DxC 700 AU)

REF B08179

For *in vitro* diagnostic use by laboratory professionals.

This document describes the general use and the instrument specific settings of the product above.

Intended purpose

The Gentian Cystatin C Immunoassay is an immunoturbidimetric assay intended for the *in vitro* quantitative determination of cystatin C in human serum and plasma on automated clinical analysers by laboratory professional users. The measurement of cystatin C is used in the diagnosis and treatment of renal diseases.

Summary and explanation of test

The non-glycosylated basic protein, cystatin C (molecular weight 13.2 kD), is produced at a constant rate in nearly every nucleated cell in the human body [1]. It is freely filtered through a normal glomerular membrane and is then reabsorbed and almost entirely catabolized in the proximal tubules. Hence, the cystatin C concentration in human blood is closely related to Glomerular Filtration Rate (GFR) [2]. A reduction in the GFR causes a rise in the concentration of cystatin C. The cystatin C concentration has not been shown to be significantly influenced by other factors such as muscular mass, inflammatory diseases, gender, age or diet [2, 3, 4].

Calibrator standardisation

The Gentian Cystatin C Calibrator is standardised against the international calibrator standard ERM-DA471/IFCC.

Relevant calculations

GFR prediction calculation

Several cystatin C based prediction equations for calculation of GFR for adults and children have been published. It should be noted that these formulas were evaluated with different cystatin C assays (Particle-Enhanced Nephelometric Immunoassay PENIA or Particle-Enhanced Turbidimetric Immunoassay PETIA) and may reveal inaccurate GFR results if an inappropriate combination of formula and assay is used. For calculation of GFR from cystatin C values measured with the Gentian assay the following prediction equation is recommended using mg/L as the unit factor [5]. The equation is valid for persons above 14 years.

$$\text{GFR [mL/min/1.73 m}^2\text{]} = \frac{79.901}{\text{Cystatin C (mg/L)}^{1.4389}}$$

Assay principle

The Gentian Cystatin C Immunoassay is a Particle-Enhanced Turbidimetric Immunoassay (PETIA). The plasma or serum sample is mixed with cystatin C immunoparticles. Cystatin C from the sample and the anti-cystatin c antibodies from the immunoparticle solution bind to form aggregates that increase the turbidity of the solution. The degree of turbidity is proportional to the concentration of cystatin C, which can be quantified via an established standard calibration curve.

Assay kit components

Products provided	BCI REF	Gentian REF
Gentian Cystatin C Reagent Kit for Beckman Coulter® AU Systems <ul style="list-style-type: none">R1 Assay Buffer (58 mL)R2 Immunoparticles (10 mL)	B08179	1103
Products required but not provided		
Gentian Cystatin C Calibrator Kit (6 levels x 1 mL)	A52763	1051
Gentian Cystatin C Control Kit (2 levels x 1 mL)	A52765	1019

All products are ready for use.

Composition

Reaction Buffer 1 (R1, 58 mL inactive ingredient): Gentian Cystatin C Assay Buffer. R1 is a MOPS [3-(N-Morpholino)-propane sulfonic acid] buffered saline, containing avian proteins and preserved with sodium azides (0.09 % (w/v)).

Reaction Buffer 2 (R2, 10 mL active ingredient): Gentian Cystatin C Immunoparticles. R2 contains a purified immunoglobulin fraction directed against human cystatin C, which is covalently attached to polystyrene nanoparticles. The solution is preserved with 0.09 % (w/v) sodium azide and antibiotics.

Warnings and precautions

1. Contains substances from human or animal origin and should be considered as potentially infectious material. Handle with caution and discard following local regulations.
2. The sodium azide concentration of the assay is not characterised as hazardous. Although, accumulated NaN_3 in lead and copper pipes may cause generation of explosive metal azides. To prevent this, rinse thoroughly if discarded into the drain.
3. Contains a sensitizing substance below concentration limit. May produce an allergic reaction in some people and may cause respiratory irritation if inhaled.
4. Contains antibiotics and must be handled with due caution.
5. Exposure may result in irritation of skin and eyes.
6. Avoid contact with incompatible materials.
7. Avoid exposure to heat and direct sunlight.

For additional safety information, please refer to the SDS (Safety Data Sheet) available on www.gentian.com.

Additional handling instructions

1. This test is for *in vitro* use only and must be handled by laboratory professionals.
2. Use only validated and approved instrument applications.
3. Do not use products after the expiration date has passed.
4. Do not mix reagents of different lots or interchange caps of reagents, controls, calibrators, and lots.
5. Tighten caps carefully back on after use of reagents, calibrators, and controls to avoid evaporation.

Cystatin C

Reagent storage and stability

All products provided for the Gentian Cystatin C Immunoassay must be stored at 2-8 °C. The expiry date is printed on the labels. The in-use stability of the Gentian Cystatin C Reagent Kit was found to be at least 9 weeks on an AU400 instrument performed as an on board study.

Specimen collection and handling

Required sample material is human serum or plasma. It is recommended to analyse the samples as fresh as possible. Sample stability testing showed that cystatin C in serum and plasma samples are stable for 14 days at room temperature (8-25 °C) and 21 days if stored at 2-8 °C. If stored below -70 °C the samples are stable for at least 5 years [6]. Mix samples well before analysing. The samples can be shipped without special cooling and must then be analysed within 14 days after shipment.

Performance characteristics

Performance characteristics AU5800

All results refer to validation of the Gentian Cystatin C Immunoassay on an AU5800 instrument at one site with one lot of reagents, unless otherwise stated.

Measuring range

The measuring range of the Gentian Cystatin C Immunoassay was found to be 0.49-7.07 mg/L. The exact measuring range is specific to the calibrator, please refer to the analytical value sheet for the lot specific calibrator values available on www.gentian.com.

Analytical sensitivity

The analytical sensitivity of the Gentian Cystatin C Immunoassay was tested based on the CLSI guideline EP17 [7]. The Limit of Quantification (LoQ) is defined as the lowest concentration of an analyte that can be reliably detected and at which the total error meets the requirements for accuracy. The LoQ of the Gentian Cystatin C Immunoassay was found to be 0.23 mg/L.

Linearity

The linearity range of the Gentian Cystatin C Immunoassay was found to be 0.49–7.07 mg/L in a linearity study based on the CLSI guideline EP06 [8].

Security zone

No antigen excess effect in samples below 32 mg/L was observed for the Gentian Cystatin C Immunoassay in a study based on the CLSI guideline EP34 [9]. Samples with a cystatin C concentration above the highest calibrator and up to 32 mg/L return a value above the highest calibrator and are flagged for rerun with automatic dilution.

Precision

Precision of the Gentian Cystatin C Immunoassay was tested in a study based on the CLSI guideline EP05 [10]. 3 serum pools and 2 controls were measured 2 times with 2 replicates (n=20).

Sample ID	Mean [mg/L]	Within run CV [%]	Between run CV [%]	Total CV [%]
P1	0.90	0.82	1.78	1.96
P2	5.29	0.49	2.05	2.10
P3	2.08	0.43	1.56	1.62
CL	0.86	1.10	3.24	3.42
CH	2.91	0.81	2.26	2.40

Recovery

Recovery was analysed by spiking a low analyte sample with a high analyte sample according to Westgard [11]. The Gentian Cystatin C Immunoassay had a recovery of 96-100 %.

Analytical specificity and limitations

There is no interference detected with the drugs tested on the recommendations from Sonntag and Scholer [12]. As the antibodies in the Gentian Cystatin C Immunoassay are of avian origin, there is no interference due to Rheumatoid Factor in the samples [13]. Interference was tested in a study based on the CLSI guideline EP07[14]. No clinically relevant difference was detected at the tested interferent concentrations.

Potential interferents	Concentration with no interference
Haemoglobin	6 g/L
Intralipid	10 g/L
Bilirubin	0.4 g/L

Instrument variation

Results obtained with the Gentian Cystatin C Immunoassay on the AU5800 instrument were compared using Passing-Bablok regression with results from AU400 instrument (Beckman Coulter) in a study based on the CLSI guideline EP09 [15].

n	Range of samples [mg/L]	Term	Co-efficient	95% CI
32	0.75-4.06	Intercept	-0.05	[-0.08, -0.02]
		Slope	1.02	[1.00, 1.06]

Performance characteristics AU680

All results refer to validation of the Gentian Cystatin C Immunoassay on an AU680 instrument at one site with one lot of reagents, unless otherwise stated.

Measuring range

The measuring range of the Gentian Cystatin C Immunoassay was found to be 0.44–7.30 mg/L. The exact measuring range is specific to the calibrator, please refer to the analytical value sheet for the lot specific calibrator values available on www.gentian.com.

Analytical sensitivity

The analytical sensitivity of the Gentian Cystatin C Immunoassay was tested in a study based on the CLSI guideline EP17 [7]. The Limit of Quantification (LoQ) is defined as the lowest concentration of an analyte that can be reliably detected and at which the total error meets the requirements for accuracy. The LoQ of the Gentian Cystatin C Immunoassay was found to be 0.28 mg/L.

Linearity

The linearity range of the Gentian Cystatin C Immunoassay was found to be 0.44–7.30 mg/L in a linearity study based on the CLSI guideline EP06 [8].

Security zone

No antigen excess effect in samples below 12 mg/L was observed for the Gentian Cystatin C Immunoassay in a study based on the CLSI guideline EP34 [9]. Samples with a cystatin C concentration above the highest calibrator and up to 12 mg/L return a value above the highest calibrator and are flagged for rerun with automatic dilution.

Precision

Precision of the Gentian Cystatin C Immunoassay was tested in study based on the CLSI guideline EP5 [10]. 4 serum pools and 2 controls were measured 2 times with 2 replicates (n=20).

Cystatin C

Sample ID	Mean [mg/L]	Within run CV [%]	Between run CV [%]	Total CV [%]
P1	0.75	0.79	2.08	2.44
P2	1.96	0.43	1.73	1.88
P3	0.80	1.09	1.35	2.00
P4	4.98	0.67	1.00	1.57
CL	1.07	0.42	1.66	2.26
CH	3.28	0.25	1.00	1.51

Recovery

Recovery was analysed by spiking a low analyte sample with a high analyte sample according to Westgard [11]. The Gentian Cystatin C Immunoassay had a recovery of 86-92 %.

Analytical specificity and limitations

There is no interference detected with the drugs tested on the recommendations from Sonntag and Scholer [12]. As the antibodies in the Gentian Cystatin C Immunoassay are of avian origin, there is no interference due to Rheumatoid Factor in the samples [13]. Interference was tested in a study based on the CLSI guideline EP07 [14]. No clinically relevant difference was detected at the tested interferent concentrations.

Potential interferents	Concentration with no interference
Haemoglobin	8.5 g/L
Intralipid	16 g/L
Bilirubin	0.2 g/L

Instrument variation

Results obtained with the Gentian Cystatin C Immunoassay on the AU680 instrument were compared using Passing-Bablok regression with results from the AU400 instrument (Beckman Coulter) in a study based on the CLSI guideline EP09 [15].

n	Range of samples [mg/L]	Term	Co-efficient	95% CI
32	0.79-4.83	Intercept	-0.02	[-0.04, 0.07]
		Slope	1.03	[0.96, 1.05]

Performance characteristics AU480

All results refer to validation of the Gentian Cystatin C Immunoassay on an AU480 instrument at one site with one lot of reagents, unless otherwise stated.

Measuring Range

The measuring range of the Gentian Cystatin C Immunoassay was found to be 0.43–7.32 mg/L. The exact measuring range is specific to the calibrator, please refer to the analytical value sheet for the lot specific calibrator values available on www.gentian.com.

Analytical sensitivity

The analytical sensitivity of the Gentian Cystatin C Immunoassay was tested in a study based on the CLSI guideline EP17 [7]. The Limit of Quantification (LoQ) is defined as the lowest concentration of an analyte that can be reliably detected and at which the total error meets the requirements for accuracy. The LoQ of the Gentian Cystatin C Immunoassay was found to be 0.43 mg/L.

Linearity

The linearity range of the Gentian Cystatin C Immunoassay was found to be 0.40–7.32 mg/L in a linearity study based on the CLSI guideline EP06 [8].

Security zone

No antigen excess effect in samples below 9.4 mg/L was observed for the Gentian Cystatin C Immunoassay in a study based on the CLSI guideline EP34 [9]. Samples with a cystatin C concentration above the highest calibrator and up to 9.4 mg/L return a value above the highest calibrator and are flagged for rerun with automatic dilution.

Precision

Precision of the Gentian Cystatin C Immunoassay was tested in study based on the CLSI guideline EP5 [10]. 3 serum pools and 2 controls were measured 2 times with 2 replicates (n=12).

Sample ID	Mean [mg/L]	Within run CV [%]	Between run CV [%]	Total CV [%]
P1	1.09	1.57	1.21	3.60
P2	3.65	0.67	0.62	1.82
P3	1.24	1.73	0.00	3.47
CL	0.87	3.10	0.00	3.72
CH	3.39	1.18	0.94	3.03

Recovery

Recovery was analysed by spiking a low analyte sample with a high analyte sample according to Westgard [11]. The Gentian Cystatin C Immunoassay had a recovery of 90-95 %.

Analytical specificity and limitations

There is no interference detected with the drugs tested on the recommendations from Sonntag and Scholer [12]. As the antibodies in the Gentian Cystatin C Immunoassay are of avian origin, there is no interference due to Rheumatoid Factor in the samples [13]. Interference was tested in a study based on the CLSI guideline EP07 [14]. No clinically relevant difference was detected at the tested interferent concentrations.

Potential interferents	Concentration with no interference
Haemoglobin	10 g/L
Intralipid	15 g/L
Bilirubin	0.6 g/L

Instrument variation

Results obtained with the Gentian Cystatin C Immunoassay on the AU480 instrument were compared using Passing-Bablok regression with results from the Architect c16000 instrument (Abbott Laboratories) in a study based on the CLSI guideline EP09 [15].

n	Range of samples [mg/L]	Term	Coefficient	95% CI
40	0.71-6.38	Intercept	0.03	[0.01, 0.04]
		Slope	0.95	[0.94, 0.97]

Cystatin C

Performance characteristics DxC 500 AU

All results refer to validation of the Gentian Cystatin C Immunoassay on a DxC 500 AU instrument at one site with one lot of reagents, unless otherwise stated.

Measuring range

The measuring range of the Gentian Cystatin C Immunoassay was found to be 0.38–7.84 mg/L. The exact measuring range is specific to the calibrator, please refer to the analytical value sheet for the lot specific calibrator values available on www.gentian.com.

Analytical sensitivity

The analytical sensitivity of the Gentian Cystatin C Immunoassay was tested in a study based on the CLSI guideline EP17 [7]. The Limit of Quantification (LoQ) is defined as the lowest concentration of an analyte that can be reliably detected and at which the total error meets the requirements for accuracy. The LoQ of the Gentian Cystatin C Immunoassay was found to be 0.32 mg/L.

Linearity

The linearity range of the Gentian Cystatin C Immunoassay was found to be 0.38–7.84 mg/L in a linearity study based on the CLSI guideline EP06 [8].

Security zone

No antigen excess effect in samples below 25.7 mg/L was observed for the Gentian Cystatin C Immunoassay in a study based on the CLSI guideline EP34 [9]. Samples with a cystatin C concentration above the highest calibrator and up to 25.7 mg/L return a value above the highest calibrator and are flagged for rerun with automatic dilution.

Precision

Precision of the Gentian Cystatin C Immunoassay was tested in a 20-day precision study based on the CLSI guideline EP05 [10]. 3 serum pools and 2 controls were measured 2 times with 2 replicates (n=80).

Sample ID	Mean [mg/L]	Within run CV [%]	Between run CV [%]	Total CV [%]
P1	0.87	0.56	1.46	2.41
P2	1.60	0.80	1.63	2.43
P3	6.37	0.73	1.63	3.66
CL	1.00	0.68	0.61	2.00
CH	3.48	0.46	0.55	1.57

Recovery

Recovery was analysed by spiking a low analyte sample with a high analyte sample according to Westgard [11]. The Gentian Cystatin C Immunoassay had a recovery of 102–109 %.

Analytical specificity and limitations

There is no interference detected with the drugs tested on the recommendations from Sonntag and Scholer [12]. As the antibodies in the Gentian Cystatin C Immunoassay are of avian origin, there is no interference due to Rheumatoid Factor in the samples [13]. Interference was tested in a study based on the CLSI guideline EP07 [14]. No clinically relevant difference was detected at the tested interferent concentrations.

Potential interferents	Concentration with no interference
Haemoglobin	8 g/L
Intralipid	10 g/L
Bilirubin	0.2 g/L

Instrument variation

Results obtained with the Gentian Cystatin C Immunoassay on the DxC 500 AU instrument were compared using Passing-Bablok regression with results from the AU5800 instrument in a study based on the CLSI guideline EP09 [15].

n	Range of samples [mg/L]	Term	Coefficient	95% CI
42	0.57– 5.72	Intercept	-0.01	[-0.05, 0.03]
		Slope	1.00	[0.97, 1.04]

Performance characteristics DxC 700 AU

All results refer to validation of the Gentian Cystatin C Immunoassay on a DxC 700 AU instrument at one site with one lot of reagents, unless otherwise stated.

Measuring range

The measuring range of the Gentian Cystatin C Immunoassay was found to be 0.40–8.07 mg/L. The exact measuring range is specific to the calibrator, please refer to the analytical value sheet for the lot specific calibrator values available on www.gentian.com.

Analytical sensitivity

The analytical sensitivity of the Gentian Cystatin C Immunoassay was tested in a study based on the CLSI guideline EP17 [7]. The Limit of Quantification (LoQ) is defined as the lowest concentration of an analyte that can be reliably detected and at which the total error meets the requirements for accuracy. The LoQ of the Gentian Cystatin C Immunoassay was found to be 0.40 mg/L.

Linearity

The linearity range of the Gentian Cystatin C Immunoassay was found to be 0.40–8.07 mg/L in a linearity study based on the CLSI guideline EP06 [8].

Security zone

No antigen excess effect in samples below 10 mg/L was observed for the Gentian Cystatin C Immunoassay in a study based on the CLSI guideline EP34 [9]. Samples with a cystatin C concentration above the highest calibrator and up to 10 mg/L return a value above the highest calibrator and are flagged for rerun with automatic dilution.

Precision

Precision of the Gentian Cystatin C Immunoassay was tested in a study based on the CLSI guideline EP5 [10]. 3 serum pools and 2 controls were measured 2 times with 2 replicates (n=80).

Sample ID	Mean [mg/L]	Within run CV [%]	Between run CV [%]	Total CV [%]
P1	0.73	0.58	0.00	0.75
P2	1.70	0.49	0.28	0.59
P3	6.13	0.44	0.18	0.60
CL	0.91	0.67	0.60	1.04
CH	3.44	0.39	0.81	0.90

Recovery

Recovery was analysed by spiking a low analyte sample with a high analyte sample according to Westgard [11]. The Gentian Cystatin C Immunoassay had a recovery of 104–105 %.

Analytical specificity and limitations

There is no interference detected with the drugs tested on the recommendations from Sonntag and Scholer [12]. As the antibodies in the Gentian Cystatin C Immunoassay are of avian origin, there is no interference due to Rheumatoid Factor in the samples [13]. Interference

Cystatin C



was tested in a study based on the CLSI guideline EP07 [14]. No clinically relevant difference was detected at the tested interferent concentrations.

Potential interferents	Concentration with no interference
Haemoglobin	10 g/L
Intralipid	10 g/L
Bilirubin	0.2 g/L

Instrument variation

Results obtained with the Gentian Cystatin C Immunoassay on the DxC 700 AU instrument were compared using Passing-Bablok regression with results from the AU5800 instrument and the Architect c4000 instrument (Abbott Laboratories) in a study based on the CLSI guideline EP09 [15].

Instru-ment	n	Range of samples [mg/L]	Term	Co-efficient	95% CI
Architect	40	0.60-6.27	Intercept	0.02	[0.00, 0.02]
			Slope	0.96	[0.95, 0.97]
AU 5800	40	0.59-6.22	Intercept	0.00	[0.00, 0.01]
			Slope	1.00	[0.99, 1.00]

Assay procedure

A detailed instrument parameter list is available in the section "Instrument Settings" below. Instrument set up, maintenance, operation and precautions must be handled in accordance with the Beckman Coulter® AU systems' instrument manuals.

Reagent preparation

The reagents are ready for use. Mix the reagents gently before placing them into the assigned reagent positions. The reagent bottles fit directly into the instrument.

Establishment of the calibration curve

Please refer to the instruction for use of the Gentian Cystatin C Calibrator Kit REF A52763 available on www.gentian.com.

QC controls

Please refer to the instruction for use of the Gentian Cystatin C Control Kit REF A52765 available on www.gentian.com.

Measuring patient samples

When a valid calibration curve has been established and the control values are within the valid range, the plasma or serum sample may be measured. Ensure that the minimum sample volume is present in the sample cups/tubes and assay the samples according to the instructions given in the instrument manual.

Results

The results are calculated automatically by the instrument for all applications established for the Gentian Cystatin C Immunoassay. The results are presented in mg/L.

Clinical performance

Sensitivity and specificity

With an eGFR cut off value of 60 mL/min/1.73 m² cystatin C has a sensitivity of 0.94 (95 % CI: 0.90-0.96) and specificity of 0.86 (95% CI: 0.78-0.91) [16].

Reference intervals

The cystatin C reference intervals were determined in a study based on the CLSI guideline C28 [17] on an Architect ci8200 instrument (Abbott Laboratories). The reference interval was determined from a population of ostensibly healthy subjects with no history of CKD. A total of 136 samples from individuals ranging in age from 20 to 84 years were measured. The samples used were serum samples. The reference

interval was calculated non-parametrically and was determined to be 0.51-1.05 mg/L. This represents the central 95 % of the population. It is recommended that every laboratory should determine a local reference interval since values may vary depending on the population tested. In a separate study involving 850 healthy children (46 % boys, 54 % girls) in the age from 5 to 15 years, the reference range 0.51-1.05 mg/L was confirmed in all ages down to 5 years of age [18].

Additional information

For more detailed information on AU Systems, refer to the appropriate system manual. Since Beckman Coulter® does not manufacture the reagent or perform quality control or other tests on individual lots, Beckman Coulter® cannot be held responsible for the quality of the data obtained which is caused by performance of the reagent, any variation between lots of reagent, or protocol changes by the manufacturer.

Shipping damage

Please notify your distributor if this product is received damaged. For technical assistance please contact your local distributor.

Symbols key

	Temperature limit
	Use by date
	Consult Instructions for Use
	Manufacturer
	CE mark with Notified Body number
	UKCA mark
	Swiss authorized representative
	<i>In Vitro</i> Diagnostic medical device
	Lot number
	Catalogue number
	Unique Device Identifier
	Contents
	R1 Assay Buffer
	R2 Immunoparticles

Cystatin C



Bjornasveien 5
N-1596 Moss
Norway
TEL: +47 99 33 99 05
www.gentian.com



Representatives

UK Responsible Person
Emergo Consulting (UK) Limited
c/o Cr360 – UL International
Compass House, Vision Park Histon
Cambridge CB24 9BZ
United Kingdom



MedEnvoy Switzerland
Gotthardstrasse 28
6302 Zug
Switzerland

References

1. Abrahamson M et al: Biochem J 1990;268:287-94.
2. Laterza OF et al: Clin Chem 2002;48:63-99.
3. Grubb AO. Adv Clin Chem 2000;35:63-99.
4. Filler G et al: Clin Biochem 2005 ;38 :1-8.
5. Flodin M et al: Scand J Clin Lab Invest 2007;67:560-567.
6. Shlipak M.G, et al: Clinical Chemistry 57: 737-745, 2011.
7. CLSI. *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline – Second Edition*. CLSI document EP17-A2. Wayne, PA: Clinical and Laboratory Standards Institute ;2012.
8. CLSI. *Evaluation of the Linearity of Quantitative Measurement Procedures*. 2nd ed. CLSI guideline EP06. Clinical Laboratory Standards Institut;2020.
9. CLSI. *Establishing and Verifying and Extended Measuring Interval Through Specimen Dilution and Spiking*; 1st Edition. CLSI guideline EP34. Wayne, PA: Clinical and Laboratory Institute; 2018.
10. CLSI. *Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline – Third Edition*. CLSI document EP05-A3. Wayne, PA: Clinical Laboratory Standards Institute; 2014.
11. Westgard JO. Basic Method Validation, 3rd Edition. 2008; ISBN13: 9781886958258.
12. Sonntag O, Scholer A. Ann Clin Biochem 2001;38:376-85.
13. Larsson A et al: J Immunol Methods. 1988 Apr 6;108(1-2):205-8.
14. CLSI. *Interference Testing in Clinical Chemistry*. 3rd ed. CLSI guideline EP07. Wayne, PA: Clinical Laboratory Standards Institute; 2018.
15. CLSI. *Measurement Procedure Comparison and Bias Estimation Using Patient Samples*. 3rd ed. CLSI guideline EP09c. Wayne, PA: Clinical Laboratory Standards Institute; 2018.
16. Qiu X et al.: Oncotarget. 2017;8(42):72985-72999.
17. CLSI. *Defining, Establishing, and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline – Third Edition*. CLSI document C28-A3c. Wayne, PA: Clinical Laboratory Standards Institute; 2008.
18. Nitsch D, et al. Am J Kidney Dis. Jun 2011;57(6):863-72.

Serious incidents

Please notify the distributor and your competent authority if any serious incidents have occurred in relation to the device.

Modifications from previous version

- Added chapter and instrument settings for the instrument DxC 500 AU.
- Included the information about the SDS available on the Gentian website.

Date of issue

2023-10-12

For other languages visit:

www.gentian.com/products/ifu/cystatin-c/beckmancoulter

Instrument settings for the Gentian Cystatin C Immunoassay

Cystatin C AU5800 application settings

System Reagent: B08179 Reagent ID: 228

Parameters		Specific Test Parameters			
General	LIH	ISE	Calculated Test		Range
Test Name: <input type="text" value="CysC"/> < > Type: <input type="text" value="Serum"/> Operation: <input type="text" value="Yes"/>					
Sample Volume	<input type="text" value="2"/> μL	Dilution	<input type="text" value="0"/> μL	OD Limit	
Pre-Dilution Rate	<input type="text" value="1"/>	Diluent Bottle	<input type="text" value="Outside"/>	Min.OD	<input type="text"/>
Rgt. Volume R1(R1-1)	<input type="text" value="150"/> μL	Dilution	<input type="text" value="0"/> μL	Reagent OD Limit	
R1-2		Dilution	<input type="text"/>	First Low	<input type="text" value="-2.0"/> High <input type="text" value="2.0"/>
R2(R2-1)	<input type="text" value="30"/> μL	Dilution	<input type="text" value="10"/> μL	Last Low	<input type="text"/>
Common Rgt. Type		Name		Dynamic Range Low	<input type="text" value="0.49"/> High <input type="text" value="7.07"/>
Wavelength	<input type="text" value="540"/> nm	Sec.	<input type="text"/>	Correlation Factor A	<input type="text" value="1.00"/> B <input type="text" value="0.00"/>
Method	<input type="text" value="End Point"/>			Factor for Maker A	<input type="text"/>
Reaction Slope	<input type="text" value="+"/>	Onboard Stability Period	<input type="text" value="60**"/> Day		<input type="text"/>
Measuring Point1 First	<input type="text" value="13"/>	Last	<input type="text" value="27"/>	LIH Influence Check	<input type="text"/>
Measuring Point2 First		Last	<input type="text"/>	Lipemia	<input type="text"/>
Linearity Limit	<input type="text"/>			Icterus	<input type="text"/>
Lag Time Check	<input type="text"/>			Hemolysis	<input type="text"/>

Parameters		Specific Test Parameters			
General	LIH	ISE	Calculated Test		Range
Test Name: <input type="text" value="CysC"/> < > Type: <input type="text" value="Serum"/>					
Value/Flag:	<input type="text" value="#"/>				
Level		Low	<input type="text" value="#"/>	High	<input type="text" value="#"/>
Specific Ranges:	From	To	Low	High	
<input type="checkbox"/> 1. Sex	Year	Month	Year	Month	
<input type="checkbox"/> 2. #	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	
<input type="checkbox"/> 3. #	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	
<input type="checkbox"/> 4. #	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	
<input type="checkbox"/> 5. #	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	
<input type="checkbox"/> 6. #	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	
<input type="checkbox"/> 7. Standard demographics			<input type="text" value="#"/>	<input type="text" value="#"/>	
<input type="checkbox"/> 8. Not within expected values			<input type="text" value="#"/>	<input type="text" value="#"/>	
Panic Value	Low	<input type="text"/>	High	<input type="text"/>	Unit <input type="text" value="mg/L"/> Decimal Places <input type="text" value="#"/>

Parameters		Calibration Parameters			
Calibrators	Calibration Specific	STAT Table Calibration			
General	ISE				
Test Name: <input type="text" value="CysC"/> < > Type: <input type="text" value="Serum"/> Cuvette: <input type="text"/>					
<input type="checkbox"/> Use Serum Cal.					
Calibration Type:	<input type="text" value="6AB"/>	Formula:	<input type="text" value="Spline"/>	Counts:	<input type="text" value="#"/>
<Calibrator Parameters>					
Point 1:	<input type="text" value="1"/>	OD	<input type="text" value="*"/>	Factor Range	
Point 2:	<input type="text" value="2"/>		<input type="text" value="*"/>	Low	<input type="text"/>
Point 3:	<input type="text" value="3"/>		<input type="text" value="*"/>	High	<input type="text"/>
Point 4:	<input type="text" value="4"/>		<input type="text" value="*"/>		
Point 5:	<input type="text" value="5"/>		<input type="text" value="*"/>		
Point 6:	<input type="text" value="6"/>		<input type="text" value="*"/>		
Point 7:	<input type="text"/>		<input type="text"/>		
Point 8:	<input type="text"/>		<input type="text"/>		
Point 9:	<input type="text"/>		<input type="text"/>		
Point 10:	<input type="text"/>		<input type="text"/>		
<Point Cal. For Master Curves>	No. of Correction Points	<input type="text"/>	Use Master Curve	<input type="text"/>	<input type="checkbox"/> Lot Calibration
Point-1	Calibrator	OD	Conc	Low	High
Point-2	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
Stability Reagent Blank <input type="text" value="28"/> Day <input type="text" value="0"/> Hour					
Calibration <input type="text" value="28"/> Day <input type="text" value="0"/> Hour					
MB Type Factor:	<input type="text"/>	1-Point Calibration Point	<input type="text"/>	<input type="checkbox"/> with Conc-0	

User defined
 * Lot specific, see analytical value sheet available on www.gentian.com
 ** Based on results from instrument AU400 (Beckman Coulter®)

Cystatin C AU680 application settings

System Reagent: B08179 Reagent ID: 228

Parameters		Specific Test Parameters			
General	LIH	ISE	Calculated Test	Range	
Test Name: <input type="text" value="CysC"/> < > Type: <input type="text" value="Serum"/> Operation <input type="text" value="Yes"/>					
Sample Volume	<input type="text" value="2"/> μL	Dilution	<input type="text" value="0"/> μL	OD Limit	
Pre-Dilution Rate	<input type="text" value="1"/>			Min.OD	<input type="text"/>
Rgt. Volume	R1(R1-1) <input type="text" value="150"/> μL	Dilution	<input type="text" value="0"/> μL	Reagent OD Limit	
				First Low	<input type="text" value="-2.0"/>
				Last Low	<input type="text"/>
				High	<input type="text" value="2.0"/>
				High	<input type="text"/>
R2(R2-1)	<input type="text" value="30"/> μL	Dilution	<input type="text" value="10"/> μL	Dynamic Range Low	<input type="text" value="0.44"/>
				High	<input type="text" value="7.30"/>
Common Rgt. Type	<input type="text"/>	Name	<input type="text"/>	Correlation Factor A	<input type="text" value="1.00"/>
Wavelength	Pri <input type="text" value="540"/> nm	Sec.	<input type="text"/>	Factor for Maker A	<input type="text" value="B"/>
Method	<input type="text" value="End Point"/>			B	<input type="text" value="0.00"/>
Reaction Slope	<input type="text" value="+"/> ∇				
Measuring Point1 First	<input type="text" value="13"/>	Last	<input type="text" value="27"/>	Onboard Stability Period	<input type="text" value="60**"/> Day <input type="text"/> Hour
Measuring Point2 First	<input type="text"/>	Last	<input type="text"/>	LIH Influence Check	<input type="text"/>
Linearity Limit	<input type="text"/>			Lipemia	<input type="text"/>
Lag Time Check	<input type="text"/>			Icterus	<input type="text"/>
				Hemolysis	<input type="text"/>

Parameters		Specific Test Parameters					
General	LIH	ISE	Calculated Test	Range			
Test Name: <input type="text" value="CysC"/> < > Type: <input type="text" value="Serum"/>							
Value/Flag:	<input type="text" value="#"/>	Low	<input type="text"/>	High	<input type="text"/>		
Level							
Specific Ranges:	From	To	Low	High	Panic Value		
	Sex	Year	Month	Year	Month	Low	High
<input type="checkbox"/> 1.	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
<input type="checkbox"/> 2.	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
<input type="checkbox"/> 3.	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
<input type="checkbox"/> 4.	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
<input type="checkbox"/> 5.	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
<input type="checkbox"/> 6.	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
<input type="checkbox"/> 7.	No demographics		<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
<input type="checkbox"/> 8.	Not within expected values		<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
Unit	<input type="text" value="mg/L"/>	Decimal Places	<input type="text" value="#"/>				

Parameters		Calibration Parameters			
Calibrators	Calibration Specific	STAT Table Calibration			
General	ISE				
Test Name: <input type="text" value="CysC"/> < > Type <input type="text" value="Serum"/> <input type="checkbox"/> Use Serum Cal.					
Calibration Type: <input type="text" value="6AB"/> Formula: <input type="text" value="Spline"/> Counts: <input type="text" value="#"/>					
<Calibrator Parameters>					
Calibrator	OD	Conc	Low	High	Slope Check
Point 1:	<input type="text" value="1"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text" value="+"/>
Point 2:	<input type="text" value="2"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	
Point 3:	<input type="text" value="3"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	
Point 4:	<input type="text" value="4"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	
Point 5:	<input type="text" value="5"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	
Point 6:	<input type="text" value="6"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	
Point 7:	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	
Point 8:	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	
Point 9:	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	
Point 10:	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	
Allowance Range Check					
<input type="checkbox"/> Reagent Blank <input type="text"/>					
<input type="checkbox"/> Calibration <input type="text"/>					
Advanced Calibration Operation <input type="text" value="#"/>					
Interval (RB/ACAL) <input type="text" value="#"/>					
<Point Cal. For No. of Correction Points <input type="text"/> Use Master Curve <input type="checkbox"/> Lot Calibration <input type="checkbox"/>					
Master Curve>					
Calibrator	OD	Conc	Low	High	Stability
Point-1	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	Reagent Blank <input type="text" value="28"/> Day <input type="text"/>
Point-2	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	Calibration <input type="text" value="28"/> Day <input type="text"/>
MB Type Factor: <input type="text"/> 1-Point Calibration Point <input type="checkbox"/> with Conc-0 <input type="checkbox"/>					

User defined

* Lot specific, see analytical value sheet available on www.gentian.com

** Based on results from instrument AU400 (Beckman Coulter®)

Cystatin C AU480 application settings

System Reagent: B08179 Reagent ID: 228

Parameters		Specific Test Parameters			
General	LIH	ISE	Calculated Test	Range	
Test Name: <input type="text" value="CysC"/> < > Type: <input type="text" value="Serum"/> Operation <input type="text" value="Yes"/>					
Sample Volume	<input type="text" value="2"/> μL	Dilution	<input type="text" value="0"/> μL	OD Limit	
Pre-Dilution Rate	<input type="text" value="1"/>			Min.OD	<input type="text"/>
Rgt. Volume	R1(R1-1) <input type="text" value="150"/> μL	Dilution	<input type="text" value="0"/> μL	Reagent OD Limit	
				First Low	<input type="text" value="-2.0"/>
				Last Low	<input type="text"/>
				High	<input type="text" value="2.0"/>
				High	<input type="text"/>
R2(R2-1)	<input type="text" value="30"/> μL	Dilution	<input type="text" value="10"/> μL	Dynamic Range Low	<input type="text" value="0.43"/>
				High	<input type="text" value="7.32"/>
Common Rgt. Type	<input type="text"/>	Name	<input type="text"/>	Correlation Factor A	<input type="text" value="1.00"/>
Wavelength	Pri <input type="text" value="540"/> nm	Sec.	<input type="text"/>	Factor for Maker A	<input type="text" value="B"/>
Method	<input type="text" value="End Point"/>			B	<input type="text" value="0.00"/>
Reaction Slope	<input type="text" value="+"/> ∇				
Measuring Point1 First	<input type="text" value="13"/>	Last	<input type="text" value="27"/>	Onboard Stability Period	<input type="text" value="60**"/> Day <input type="text"/> Hour
Measuring Point2 First	<input type="text"/>	Last	<input type="text"/>	LIH Influence Check	<input type="text"/>
Linearity Limit	<input type="text"/>			Lipemia	<input type="text"/>
Lag Time Check	<input type="text"/>			Icterus	<input type="text"/>
				Hemolysis	<input type="text"/>

Parameters		Specific Test Parameters					
General	LIH	ISE	Calculated Test	Range			
Test Name: <input type="text" value="CysC"/> < > Type: <input type="text" value="Serum"/>							
Value/Flag:	<input type="text" value="#"/>	Low	<input type="text"/>	High	<input type="text"/>		
Level							
Specific Ranges:	From	To	Low	High	Panic Value		
	Sex	Year	Month	Year	Month	Low	High
<input type="checkbox"/> 1.	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
<input type="checkbox"/> 2.	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
<input type="checkbox"/> 3.	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
<input type="checkbox"/> 4.	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
<input type="checkbox"/> 5.	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
<input type="checkbox"/> 6.	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
<input type="checkbox"/> 7.	No demographics		<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
<input type="checkbox"/> 8.	Not within expected values		<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
Unit	<input type="text" value="mg/L"/>	Decimal Places	<input type="text" value="#"/>				

Parameters		Calibration Parameters			
Calibrators	Calibration Specific	STAT Table Calibration			
General	ISE				
Test Name: <input type="text" value="CysC"/> < > Type <input type="text" value="Serum"/> <input type="checkbox"/> Use Serum Cal.					
Calibration Type: <input type="text" value="6AB"/> Formula: <input type="text" value="Spline"/> Counts: <input type="text" value="#"/>					
<Calibrator Parameters>					
Calibrator	OD	Conc	Low	High	Slope Check
Point 1:	<input type="text" value="1"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text" value="+"/> ∇
Point 2:	<input type="text" value="2"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	Allowance Range Check
Point 3:	<input type="text" value="3"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="checkbox"/> Reagent Blank <input type="text"/>
Point 4:	<input type="text" value="4"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="checkbox"/> Calibration <input type="text"/>
Point 5:	<input type="text" value="5"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	Advanced Calibration
Point 6:	<input type="text" value="6"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	Operation <input type="text" value="#"/>
Point 7:	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	Interval (RB/ACAL) <input type="text" value="#"/>
Point 8:	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	
Point 9:	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	
Point 10:	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	
<Point Cal. For No. of Correction Points <input type="text"/> Use Master Curve <input type="checkbox"/> Lot Calibration <input type="checkbox"/>					
Master Curve>					
Calibrator	OD	Conc	Low	High	Stability
Point-1	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	Reagent Blank <input type="text" value="28"/> Day <input type="text"/> Hour
Point-2	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	Calibration <input type="text" value="28"/> Day <input type="text"/> Hour
MB Type Factor: <input type="text"/> 1-Point Calibration Point <input type="checkbox"/> with Conc-0 <input type="checkbox"/>					

User defined

* Lot specific, see analytical value sheet available on www.gentian.com

** Based on results from instrument AU400 (Beckman Coulter®)

Cystatin C



Cystatin C DxC 500 AU application settings

System Reagent: B08179 Reagent ID: 228

TEST CONFIGURATION & CHEMISTRY DETAILS																	
Assay Name	Test	Rev	Discipline		Chemistry												
Test ID	CYS		Calculated Result		<input type="checkbox"/>												
LIS Code	CYS																
UNITS AND RANGE SETTINGS																	
Use Settings from	Serum	Units	mg/L	Decimal Places	x.xx												
Test Kind	General	Revision	01	<input checked="" type="checkbox"/> Multi Reagent Switch													
Reagent Name	CYS	Reagent ID	228	<input type="checkbox"/> FSE Test													
ABB Name	CYS1G	Parameter Long Name															
Cystatin C B08179 CYS1G CYS1G Serum																	
Region	<input checked="" type="checkbox"/> US	<input checked="" type="checkbox"/> OUS	<input checked="" type="checkbox"/> AP	<input type="checkbox"/> JP	<input checked="" type="checkbox"/> EU	<input type="checkbox"/> Other											
GENERAL PARAMETERS																	
SAMPLE VOLUME			Sample Volume	2.0	µL	Dilution	0	µL	REACTION OD LIMIT								
REAGENT VOLUME			Predilution Rate	1				Low		High							
WAVELENGTH			R1-1	150	µL	Dilution	0	µL	REACTION BLANK OD LIMIT								
METHOD			R2-1	30	µL	Dilution	10	µL	First: Low	-2.0000	High	2.0000					
REACTION SLOPE									Last: Low	-2.0000	High	2.0000					
MEASURING POINT									ANALYTICAL MEASURING RANGE								
LINEARITY LIMIT									Low	0.38	High	7.84					
LAG TIME CHECK									MANUFACTURER FACTOR								
PERFORM LIH CHECK									A	1	B	0					
PERFORM LAG TIME CHECK									REAGENT ONBOARD STABILITY								
PERFORM LAG TIME CHECK									60** Days			0** Hours					
PERFORM LAG TIME CHECK									LIH INFLUENCE CHECK								
PERFORM LAG TIME CHECK									<input type="checkbox"/> Perform LIH check								
PERFORM LAG TIME CHECK									Lipemia			+ ▼					
PERFORM LAG TIME CHECK									Icterus			+ ▼					
PERFORM LAG TIME CHECK									Hemolysis			+ ▼					
CALIBRATION PARAMETERS																	
Base Unit	Decimal Place	Unit 1	Factor 1	Unit 2	Factor 2	Unit 3	Factor 3	Unit 4	Factor 4								
mg/L	2	None	0	None	0	None	0	None	0								
CALIBRATOR SPECIFIC					CALIBRATION OD AND CONCENTRATION PARAMETERS												
Calibration Type					Counts	<input type="checkbox"/> Use highest calibrator for Upper AMR											
Formula					MB Factor	Calibrator Name				Conc	OD Range Low	OD Range High					
Calibrator Name					Positive Cutoff	Point 1	CYS CAL-1	*	-2.0000	2.0000							
<input checked="" type="checkbox"/> SLOPE CHECK					Number of Levels	Point 2	CYS CAL-2	*	-2.0000	2.0000							
Slope Check						Point 3	CYS CAL-3	*	-2.0000	2.0000							
STABILITY AND INTERVAL						Point 4	CYS CAL-4	*	-2.0000	2.0000							
Reagent Blank Stability					28 Days	0 Hours	Interval	Lot	Point 5	CYS CAL-5	*	-2.0000	2.0000				
Calibration Stability					28 Days	0 Hours	Interval	Lot	Point 6	CYS CAL-6	*	-2.0000	2.0000				
									Point 7								
									OD DELTA CHECK								
									<input type="checkbox"/> Reagent Blank				0.0000				
									<input type="checkbox"/> Calibration				0.0000				
PROZONE CHECK PARAMETERS																	
<input type="checkbox"/> Logic Check 1				<input type="checkbox"/> Logic Check 2				<input type="checkbox"/> Logic Check 3									
Check Points			Decision Values			Check Points			Decision Values			Check Points			Decision Values		
Point 1	#		Value 1	#		Point 1	#		Value 1	#		Point 1	#		Value 1	#	
Point 2	#		Value 2	#		Interval	#		Value 2	#		Interval	#		Value 2	#	
Point 3	#		Value 3	#		Limit Points			Limit Points			Limit Points					
Limit 1	#					Limit 1	#					Limit 1	#				
Limit 2	#					Limit 2	#					Limit 2	#				
Check Pattern																	
Pattern																	

User defined

* Lot specific, see analytical value sheet available on www.gentian.com

** Based on results from instrument AU400 (Beckman Coulter®)

Cystatin C



Cystatin C DxC 700 AU application settings

System Reagent: B08179 Reagent ID: 228

General	LIH	ISE	Calculated Test	Range
Test Name: CYS1G Test No Type: Serum Operation: Yes				
Sample Volume	2.0 µL	Dilution 0 µL	OD Limit	
Pre-Dilution Rate	1		Reagent OD Limit	Min. OD Max OD
Reagent Volume R1 (R1-1)	150 µL	Dilution 0 µL	1st	Low -2.0000 High 2.0000
R1-2		Dilution	Last	Low -2.0000 High 2.0000
R2 (R2-1)	30 µL	Dilution 10 µL	Analytical Measuring Range	Low 0.40 High 8.07
Common Reagent Type	None	Name None	Correlation Factor	A 1 B 0
Wavelength Pri	540 nm	Sec None nm	Manufacturer Factor	A 1 B 0
Method	END			
Reaction Slope	+		Onboard Stability Period	60** Day 0 Hour
Measuring Point-1 1st	13	Last 27	LIH Influence Check	No
Measuring Point-2 1st		Last	Lipemia	
Linearity Limit			Icterus	
Lag Time Check			Hemolysis	

General	LIH	ISE	Calculated Test	Range					
Test Name: CYS1G Test No Type: Serum									
Value/Flag	Value	Level	Low -9999.99	High 9999.99					
Specific Ranges									
	Sex	Year	Month	Year	To	Month	Other Type	Low	High
<input type="checkbox"/> 1:	#	#	#	#	#	#	None	#	#
<input type="checkbox"/> 2:	#	#	#	#	#	#	None	#	#
<input type="checkbox"/> 3:	#	#	#	#	#	#	None	#	#
<input type="checkbox"/> 4:	#	#	#	#	#	#	None	#	#
<input type="checkbox"/> 5:	#	#	#	#	#	#	None	#	#
<input type="checkbox"/> 6:	#	#	#	#	#	#	None	#	#
7:	Standard demographics							#	#
8:	Not within expected values							#	#
Critical Limits	Low #	High #	Unit mg/L	Select	Decimal Places 2				

Calibrators	General	ISE			
Test Name: CYS1G Type: Serum					
<input type="checkbox"/> Use Serum Cal.					
Calibration Type:	6AB	Formula: Spline			
Counts:	2				
<Calibrator Parameters>					
	Calibrator	OD	Conc	Range	
				Low	High
Point-1	CYSC Calibrator Level 1		*	-2.0000	2.0000
Point-2	CYSC Calibrator Level 2		*	-2.0000	2.0000
Point-3	CYSC Calibrator Level 3		*	-2.0000	2.0000
Point-4	CYSC Calibrator Level 4		*	-2.0000	2.0000
Point-5	CYSC Calibrator Level 5		*	-2.0000	2.0000
Point-6	CYSC Calibrator Level 6		*	-2.0000	2.0000
Point-7					
MB Type Factor		1-Point Calibration Point	None	<input type="checkbox"/> with Conc-0	
				Stability	
				Reagent Blank	28 Day 0 Hour
				Calibration	28 Day 0 Hour
				Slope Check	+
				Allowable Range Check	
				<input type="checkbox"/> Reagent Blank	
				<input type="checkbox"/> Calibration	
				Advanced Calibration	
				Operation	No
				Interval (RB)	
				Interval (ACAL)	

User defined

* Lot specific, see analytical value sheet available on www.gentian.com

** Based on results from instrument AU400 (Beckman Coulter®)