Cystatin C

Gentian Cystatin C Immunoassay pe sistemele Beckman Coulter® AU (AU5800, AU680, AU480, DxC 500 AU/DxC 500i AU, DxC 700 AU)

REF B08179

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R2 Imuno-particule

Bjornasveien 5
N-1596 Moss
Norvegia
TEL: +47 99 33 99 05
www.gentian.com

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Emergo Consulting (UK) Limited
c/o Cr360 – UL International
Compass House, Vision Park Histon
Cambridge CB24 9BZ
Regatul Unit

MedEnvoy Switzerland
Gotthardstrasse 28
6302 Zug
Elveția

Referințe

1. Abrahamson M et al: Biochem J 1990;268:287-94.
2. Laterza OF et al: Clin Chem 2002;48:63-99.
3. Grubb AO. Adv Clin Chem 2000;35:63-99.

4. Filler G et al: Clin Biochem 2005 ;38 :1-8.
5. Flodin M et al: Scand J Clin Lab Invest 2007;67:560-567.
6. Shlipak M.G, et al: Clinical Chemistry 57: 737-745, 2011.
7. CLSI. *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline – Second Edition.* CLSI document EP17-A2. Wayne, PA: Clinical and Laboratory Standards Institute ;2012.
8. CLSI. *Evaluation of the Linearity of Quantitative Measurement Procedures.* 2nd ed. CLSI guideline EP06. Clinical Laboratory Standards Institut;2020.
9. CLSI. *Establishing and Verifying and Extended Measuring Interval Through Specimen Dilution and Spiking;* 1st Edition. CLSI guideline EP34. Wayne, PA: Clinical and Laboratory Institute; 2018.
10. CLSI. *Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline – Third Edition.* CLSI document EP05-A3. Wayne, PA: Clinical Laboratory Standards Institute; 2014.
11. Westgard JO. Basic Method Validation, 3rd Edition. 2008; ISBN13: 9781886958258.
12. Sonntag O, Scholer A. Ann Clin Biochem 2001;38:376-85.
13. Larsson A et al: J Immunol Methods. 1988 Apr 6;108(1-2):205-8.
14. CLSI. *Interference Testing in Clinical Chemistry.* 3rd ed. CLSI guideline EP07. Wayne, PA: Clinical Laboratory Standards Institute; 2018.
15. CLSI. *Measurement Procedure Comparison and Bias Estimation Using Patient Samples.* 3rd ed. CLSI guideline EP09c. Wayne, PA: Clinical Laboratory Standards Institute; 2018.
16. Qiu X et al.: Oncotarget. 2017;8(42):72985-72999.
17. CLSI. *Defining, Establishing, and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline – Third Edition.* CLSI document C28-A3c. Wayne, PA: Clinical Laboratory Standards Institute; 2008.
18. Nitsch D, et al. Am J Kidney Dis. Jun 2011;57(6):863-72.

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- Inclus DxC 500i.

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2025-04-03

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www.gentian.com/products/ifu/cystatin-c/beckmancoulter

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Reactiv de sistem: B08179 Date de identificare (ID) reactiv: 228

TEST CONFIGURATION & CHEMISTRY DETAILS									
Assay Name	Test	Rev	Discipline		Chemistry				
Test ID	CYS		Calculated Result		☐				
LIS Code	CYS								
UNITS AND RANGE SETTINGS									
Use Settings from	Serum	Units	mg/L	Decimal Places	x.xx	Plasma			
Test Kind	General	Revision	01	☒ Multi Reagent Switch					
Reagent Name	CYS	Reagent ID	228	☐ FSE Test					
ABB Name	CYS1G	Parameter Long Name	Cystatin C B08179 CYS1G CYSC Serum						
Region	☒ US	☒ OUS	☒ AP	☐ JP	☒ EU	☐ Other			
GENERAL PARAMETERS									
SAMPLE VOLUME		Sample Volume	2.0	µL	Dilution	0	µL	REACTION OD LIMIT	
REAGENT VOLUME		Predilution Rate	1				Low		High
WAVELENGTH		R1-1	150	µL	Dilution	0	µL	REACTION BLANK OD LIMIT	
METHOD		R2-1	30	µL	Dilution	10	µL	First: Low	
REACTION SLOPE								Last: Low	
MEASURING POINT								High	
Linearity Limit								ANALYTICAL MEASURING RANGE	
Lag Time Check								Low	
								High	
								MANUFACTURER FACTOR	
								A	
								B	
								REAGENT ONBOARD STABILITY	
								60** Days	
								0** Hours	
								LIH INFLUENCE CHECK	
								☐ Perform LIH check	
								Lipemia	
								Icterus	
								Hemolysis	
CALIBRATION PARAMETERS									
Base Unit	Decimal Place	Unit 1	Factor 1	Unit 2	Factor 2	Unit 3	Factor 3	Unit 4	Factor 4
mg/L	2	None	0	None	0	None	0	None	0
CALIBRATOR SPECIFIC									
Calibration Type		6AB		Counts		2		☐ Use highest calibrator for Upper AMR	
Formula		Spline		MB Factor					
Calibrator Name		Add CYS		Positive Cutoff					
☒ SLOPE CHECK		Number of Levels		6					
Slope Check		+							
STABILITY AND INTERVAL									
Reagent Blank Stability		28 Days		0 Hours		Interval		Lot	
Calibration Stability		28 Days		0 Hours		Interval		Lot	
OD DELTA CHECK									
☐ Reagent Blank		0.0000							
☐ Calibration		0.0000							
PROZONE CHECK PARAMETERS									
☐ Logic Check 1		Check Points		Decision Values		☐ Logic Check 2		Check Points	
Point 1		#		Value 1		Point 1		#	
Point 2		#		Value 2		Interval		#	
Point 3		#		Value 3					
Limit Points		Limit 1		#		Limit Points		Limit 1	
Limit 2		#				Limit 2		#	
Check Pattern		Pattern		#					

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Reactiv de sistem: B08179 Date de identificare (ID) reactiv: 228

General	LIH	ISE	Calculated Test	Range																																																																																				
Test Name: CYS1G Test No Type: Serum Operation Yes																																																																																								
<table> <tr> <td>Sample Volume</td> <td> 2.0 µL</td> <td>Dilution</td> <td> 0 µL</td> <td>OD Limit</td> <td></td> </tr> <tr> <td>Pre-Dilution Rate</td> <td> 1 </td> <td></td> <td></td> <td>Min. OD</td> <td> </td> </tr> <tr> <td>Reagent Volume</td> <td>R1 (R1-1) 150 µL</td> <td>Dilution</td> <td> 0 µL</td> <td>Reagent OD Limit</td> <td>Max OD </td> </tr> <tr> <td></td> <td>R1-2 µL</td> <td>Dilution</td> <td> µL</td> <td>1st</td> <td>Low -2.0000 High 2.0000 </td> </tr> <tr> <td></td> <td>R2 (R2-1) 30 µL</td> <td>Dilution</td> <td> 10 µL</td> <td>Last</td> <td>Low -2.0000 High 2.0000 </td> </tr> <tr> <td>Common Reagent</td> <td>Type None </td> <td>Name</td> <td> None </td> <td>Analytical Measuring Range</td> <td>Low 0.40 High 8.07 </td> </tr> <tr> <td></td> <td></td> <td></td> <td></td> <td>Correlation Factor</td> <td>A 1 B 0 </td> </tr> <tr> <td>Wavelength</td> <td>Pri 540 nm</td> <td>Sec</td> <td> None nm</td> <td>Manufacturer Factor</td> <td>A 1 B 0 </td> </tr> <tr> <td>Method</td> <td> END </td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td>Reaction Slope</td> <td> + ▼</td> <td></td> <td></td> <td>Onboard Stability Period</td> <td> 60 ** Day 0 Hour</td> </tr> <tr> <td>Measuring Point-1</td> <td>1st 13 </td> <td>Last</td> <td> 27 </td> <td>LIH Influence Check</td> <td> No ▼</td> </tr> <tr> <td>Measuring Point-2</td> <td>1st </td> <td>Last</td> <td> </td> <td>Lipemia</td> <td> ▼</td> </tr> <tr> <td>Linearity Limit</td> <td> %</td> <td></td> <td></td> <td>Icterus</td> <td> ▼</td> </tr> <tr> <td>Lag Time Check</td> <td> ▼</td> <td></td> <td></td> <td>Hemolysis</td> <td> ▼</td> </tr> </table>					Sample Volume	2.0 µL	Dilution	0 µL	OD Limit		Pre-Dilution Rate	1			Min. OD		Reagent Volume	R1 (R1-1) 150 µL	Dilution	0 µL	Reagent OD Limit	Max OD		R1-2 µL	Dilution	µL	1 st	Low -2.0000 High 2.0000		R2 (R2-1) 30 µL	Dilution	10 µL	Last	Low -2.0000 High 2.0000	Common Reagent	Type None	Name	None	Analytical Measuring Range	Low 0.40 High 8.07					Correlation Factor	A 1 B 0	Wavelength	Pri 540 nm	Sec	None nm	Manufacturer Factor	A 1 B 0	Method	END					Reaction Slope	+ ▼			Onboard Stability Period	60 ** Day 0 Hour	Measuring Point-1	1st 13	Last	27	LIH Influence Check	No ▼	Measuring Point-2	1st	Last		Lipemia	▼	Linearity Limit	%			Icterus	▼	Lag Time Check	▼			Hemolysis	▼
Sample Volume	2.0 µL	Dilution	0 µL	OD Limit																																																																																				
Pre-Dilution Rate	1			Min. OD																																																																																				
Reagent Volume	R1 (R1-1) 150 µL	Dilution	0 µL	Reagent OD Limit	Max OD																																																																																			
	R1-2 µL	Dilution	µL	1 st	Low -2.0000 High 2.0000																																																																																			
	R2 (R2-1) 30 µL	Dilution	10 µL	Last	Low -2.0000 High 2.0000																																																																																			
Common Reagent	Type None	Name	None	Analytical Measuring Range	Low 0.40 High 8.07																																																																																			
				Correlation Factor	A 1 B 0																																																																																			
Wavelength	Pri 540 nm	Sec	None nm	Manufacturer Factor	A 1 B 0																																																																																			
Method	END																																																																																							
Reaction Slope	+ ▼			Onboard Stability Period	60 ** Day 0 Hour																																																																																			
Measuring Point-1	1st 13	Last	27	LIH Influence Check	No ▼																																																																																			
Measuring Point-2	1st	Last		Lipemia	▼																																																																																			
Linearity Limit	%			Icterus	▼																																																																																			
Lag Time Check	▼			Hemolysis	▼																																																																																			

General	LIH	ISE	Calculated Test	Range																																																																																	
Test Name: CYS1G Test No Type: Serum																																																																																					
<table> <tr> <td>Value/Flag</td> <td> Value </td> <td>Level</td> <td>Low -99999.99 High 99999.99 </td> </tr> </table>					Value/Flag	Value	Level	Low -99999.99 High 99999.99																																																																													
Value/Flag	Value	Level	Low -99999.99 High 99999.99																																																																																		
Specific Ranges <table> <tr> <th></th> <th>Sex</th> <th>Year</th> <th>Month</th> <th>Year</th> <th>Month</th> <th>Other Type</th> <th>Low</th> <th>High</th> </tr> <tr> <td>□ 1:</td> <td> # </td> <td> # </td> <td> # </td> <td> # </td> <td> # </td> <td> None </td> <td> # </td> <td> # </td> </tr> <tr> <td>□ 2:</td> <td> # </td> <td> # </td> <td> # </td> <td> # </td> <td> # </td> <td> None </td> <td> # </td> <td> # </td> </tr> <tr> <td>□ 3:</td> <td> # </td> <td> # </td> <td> # </td> <td> # </td> <td> # </td> <td> None </td> <td> # </td> <td> # </td> </tr> <tr> <td>□ 4:</td> <td> # </td> <td> # </td> <td> # </td> <td> # </td> <td> # </td> <td> None </td> <td> # </td> <td> # </td> </tr> <tr> <td>□ 5:</td> <td> # </td> <td> # </td> <td> # </td> <td> # </td> <td> # </td> <td> None </td> <td> # </td> <td> # </td> </tr> <tr> <td>□ 6:</td> <td> # </td> <td> # </td> <td> # </td> <td> # </td> <td> # </td> <td> None </td> <td> # </td> <td> # </td> </tr> <tr> <td>7:</td> <td colspan="6">Standard demographics</td> <td> # </td> <td> # </td> </tr> <tr> <td>8:</td> <td colspan="6">Not within expected values</td> <td> # </td> <td> # </td> </tr> </table>						Sex	Year	Month	Year	Month	Other Type	Low	High	□ 1:	#	#	#	#	#	None	#	#	□ 2:	#	#	#	#	#	None	#	#	□ 3:	#	#	#	#	#	None	#	#	□ 4:	#	#	#	#	#	None	#	#	□ 5:	#	#	#	#	#	None	#	#	□ 6:	#	#	#	#	#	None	#	#	7:	Standard demographics						#	#	8:	Not within expected values						#	#
	Sex	Year	Month	Year	Month	Other Type	Low	High																																																																													
□ 1:	#	#	#	#	#	None	#	#																																																																													
□ 2:	#	#	#	#	#	None	#	#																																																																													
□ 3:	#	#	#	#	#	None	#	#																																																																													
□ 4:	#	#	#	#	#	None	#	#																																																																													
□ 5:	#	#	#	#	#	None	#	#																																																																													
□ 6:	#	#	#	#	#	None	#	#																																																																													
7:	Standard demographics						#	#																																																																													
8:	Not within expected values						#	#																																																																													
<table> <tr> <td>Critical Limits</td> <td>Low # </td> <td>High # </td> <td>Unit mg/L </td> <td>Select </td> </tr> <tr> <td></td> <td></td> <td></td> <td></td> <td>Decimal Places 2 </td> </tr> </table>					Critical Limits	Low #	High #	Unit mg/L	Select					Decimal Places 2																																																																							
Critical Limits	Low #	High #	Unit mg/L	Select																																																																																	
				Decimal Places 2																																																																																	

Calibrators	General	ISE																																																						
Test Name: CYS1G Type: Serum																																																								
☐ Use Serum Cal.																																																								
<table> <tr> <td>Calibration Type:</td> <td> 6AB </td> <td>Formula:</td> <td> Spline </td> <td>Counts:</td> <td> 2 </td> </tr> </table>			Calibration Type:	6AB	Formula:	Spline	Counts:	2																																																
Calibration Type:	6AB	Formula:	Spline	Counts:	2																																																			
<Calibrator Parameters> <table> <tr> <th></th> <th>Calibrator</th> <th>OD</th> <th>Conc</th> <th colspan="2">Range</th> </tr> <tr> <th></th> <th></th> <th></th> <th></th> <th>Low</th> <th>High</th> </tr> <tr> <td>Point-1</td> <td>CYSC Calibrator Level 1 ▼</td> <td></td> <td>*</td> <td>-2.0000</td> <td>2.0000</td> </tr> <tr> <td>Point-2</td> <td>CYSC Calibrator Level 2 ▼</td> <td></td> <td>*</td> <td>-2.0000</td> <td>2.0000</td> </tr> <tr> <td>Point-3</td> <td>CYSC Calibrator Level 3 ▼</td> <td></td> <td>*</td> <td>-2.0000</td> <td>2.0000</td> </tr> <tr> <td>Point-4</td> <td>CYSC Calibrator Level 4 ▼</td> <td></td> <td>*</td> <td>-2.0000</td> <td>2.0000</td> </tr> <tr> <td>Point-5</td> <td>CYSC Calibrator Level 5 ▼</td> <td></td> <td>*</td> <td>-2.0000</td> <td>2.0000</td> </tr> <tr> <td>Point-6</td> <td>CYSC Calibrator Level 6 ▼</td> <td></td> <td>*</td> <td>-2.0000</td> <td>2.0000</td> </tr> <tr> <td>Point-7</td> <td>CYSC Calibrator Level 7 ▼</td> <td></td> <td>*</td> <td>-2.0000</td> <td>2.0000</td> </tr> </table>				Calibrator	OD	Conc	Range						Low	High	Point-1	CYSC Calibrator Level 1 ▼		*	-2.0000	2.0000	Point-2	CYSC Calibrator Level 2 ▼		*	-2.0000	2.0000	Point-3	CYSC Calibrator Level 3 ▼		*	-2.0000	2.0000	Point-4	CYSC Calibrator Level 4 ▼		*	-2.0000	2.0000	Point-5	CYSC Calibrator Level 5 ▼		*	-2.0000	2.0000	Point-6	CYSC Calibrator Level 6 ▼		*	-2.0000	2.0000	Point-7	CYSC Calibrator Level 7 ▼		*	-2.0000	2.0000
	Calibrator	OD	Conc	Range																																																				
				Low	High																																																			
Point-1	CYSC Calibrator Level 1 ▼		*	-2.0000	2.0000																																																			
Point-2	CYSC Calibrator Level 2 ▼		*	-2.0000	2.0000																																																			
Point-3	CYSC Calibrator Level 3 ▼		*	-2.0000	2.0000																																																			
Point-4	CYSC Calibrator Level 4 ▼		*	-2.0000	2.0000																																																			
Point-5	CYSC Calibrator Level 5 ▼		*	-2.0000	2.0000																																																			
Point-6	CYSC Calibrator Level 6 ▼		*	-2.0000	2.0000																																																			
Point-7	CYSC Calibrator Level 7 ▼		*	-2.0000	2.0000																																																			
<table> <tr> <td>MB Type Factor</td> <td> </td> <td>1-Point Calibration Point</td> <td> None </td> <td>☐ with Conc-0</td> </tr> </table>			MB Type Factor		1-Point Calibration Point	None	☐ with Conc-0																																																	
MB Type Factor		1-Point Calibration Point	None	☐ with Conc-0																																																				
<table> <tr> <td>Slope Check</td> <td> + ▼</td> </tr> <tr> <td>Allowable Range Check</td> <td>☐ Reagent Blank </td> </tr> <tr> <td></td> <td>☐ Calibration </td> </tr> <tr> <td>Advanced Calibration</td> <td>☐ Operation No </td> </tr> <tr> <td>Interval (RB)</td> <td> </td> </tr> <tr> <td>Interval (ACAL)</td> <td> </td> </tr> <tr> <td>Stability</td> <td></td> </tr> <tr> <td>Reagent Blank</td> <td> 28 Day 0 Hour</td> </tr> <tr> <td>Calibration</td> <td> 28 Day 0 Hour</td> </tr> </table>			Slope Check	+ ▼	Allowable Range Check	☐ Reagent Blank		☐ Calibration	Advanced Calibration	☐ Operation No	Interval (RB)		Interval (ACAL)		Stability		Reagent Blank	28 Day 0 Hour	Calibration	28 Day 0 Hour																																				
Slope Check	+ ▼																																																							
Allowable Range Check	☐ Reagent Blank																																																							
	☐ Calibration																																																							
Advanced Calibration	☐ Operation No																																																							
Interval (RB)																																																								
Interval (ACAL)																																																								
Stability																																																								
Reagent Blank	28 Day 0 Hour																																																							
Calibration	28 Day 0 Hour																																																							

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