

Cystatin C



Gentian Cystatin C Immunoassay on Beckman Coulter® AU Systems (AU5800, AU680, AU480, DxC 500 AU/DxC 500i, DxC 700 AU)

REF B08179

This document describes the general use and instrument specific settings of the product above and is applicable for USA only.

Intended use

The Cystatin C Immunoassay on the Beckman Coulter® AU Systems is an *in vitro* diagnostic test for quantitative determination of cystatin C in human serum and plasma. The measurement of cystatin C is used in the diagnosis and treatment of renal diseases.

Summary and explanation of test

The non-glycosylated basic protein, cystatin C (molecular weight 13.2 kD), is produced at a constant rate in nearly every nucleated cell in the human body [1]. It is freely filtered through a normal glomerular membrane and is then reabsorbed and almost entirely catabolized in the proximal tubules. Hence, the cystatin C concentration in human blood is closely related to Glomerular Filtration Rate (GFR) [2]. A reduction in the GFR causes a rise in the concentration of cystatin C. The cystatin C concentration has not been shown to be significantly influenced by other factors such as muscular mass, inflammatory diseases, gender, age or diet [2, 3, 4].

Principles of the procedure

Serum or plasma sample from human is mixed with cystatin C immunoparticles. Cystatin C from the sample and anti-cystatin C from the immunoparticles aggregates. The complex particles created absorb light, and by turbidimetry the absorption is related to cystatin C concentration via interpolation on an established standard calibration curve. The AU platforms will automatically calculate the results.

Composition

Reaction Buffer 1 (R1, 58 mL inactive ingredient): Gentian Cystatin C Assay Buffer. R1 is a MOPS [3-(N-Morpholino)-propane sulfonic acid] buffered saline, containing avian proteins and preserved with sodium azides (0.09 % (w/v)).

Reaction Buffer 2 (R2, 10 mL active ingredient): Gentian Cystatin C Immunoparticles. R2 contains a purified immunoglobulin fraction directed against human cystatin C, which is covalently attached to polystyrene nanoparticles. The solution is preserved with 0.09 % (w/v) sodium azide and antibiotics.

Warnings and precautions

For *in vitro* diagnostic use by laboratory professionals.

Caution: Federal law restricts this device to sale by or on the order of a physician.

1. Contains substances from human or animal origin and should be considered as potentially infectious material. Serum used in Gentian Cystatin C Controls and calibrators is tested for hepatitis HBsAG, anti-HCV, anti-HIV1 and anti-HIV2 and found to be negative. Handle with caution and discard following local regulations.
2. Contains antibiotics and must be handled with due caution.
3. The sodium azide concentration of the assay is not characterized as hazardous. Although, accumulated NaN₃ in lead and copper pipes may cause generation of explosive metal azides. To prevent this, rinse thoroughly if discarded into the drain.

4. Exposure may result in irritation of skin and eyes.
5. Avoid contact with incompatible materials.
6. Avoid exposure to heat and direct sunlight.

To obtain the SDS (Safety Data Sheet), please refer to the SDS (Safety Data Sheet) available on www.gentian.com.

Additional handling instructions

1. This test is for *in vitro* use only and must be handled by laboratory professionals.
2. Use only validated and approved instrument applications.
3. Do not use products after the expiration date has passed.
4. Do not mix reagents of different lots or interchange caps of reagents, controls, calibrators, and lots.
5. Tighten caps carefully back on after use of reagents, calibrators, and controls to avoid evaporation.

Directions for reconstitution/dilution

The product is ready to use.

Storage instructions

Shelf life of unopened reagents at 2-8 °C: See expiry date on the label.

Specimen collection and preparation

The required sample material is human serum or EDTA/heparinized plasma. It is recommended to analyze the samples as fresh as possible. Sample stability testing showed that cystatin C in serum and plasma samples are stable for 14 days at room temperature (8-25 °C), 21 days if stored at 2-8 °C. If stored below -70 °C the samples are stable for at least 5 years [5]. Mix samples well before analyzing them.

Procedure

A detailed instrument parameter list is available in the section "Instrument settings for Cystatin C Immunoassay" below. Instrument set up, maintenance, operation and precautions must be handled in accordance with the Beckman Coulter® AU system's instrument manuals.

Reagent preparation

Gentian Cystatin C reagents are supplied ready for use. Mix gently before loading into instrument. Reagents should be stored capped at 2-8 °C when not in use.

Assay kit components

Materials provided	
Gentian Cystatin C Reagent Kit for Beckman Coulter® • R1 Assay Buffer (58 mL) • R2 Immunoparticles (10 mL)	REF B08179
Materials required, but not provided	
Gentian Cystatin C Calibrator Kit (6 levels x 1 mL)	REF A52763
Gentian Cystatin C Control Kit (2 levels x 1 mL)	REF A52765

Cystatin C

All materials are ready for use.

Stability

Stability after opening: Until expiry date at 2-8 °C. On-board stability: 9 weeks at correct temperature (2-8 °C).

Calibrator standardization

Gentian Cystatin C Calibrator is standardized against the international calibrator standard ERM-DA471/IFCC.

Establishment of the calibration curve

Please refer to the package insert of the Gentian Cystatin C Calibrator Kit REF A52763.

QC controls

Please refer to the package insert of the Gentian Cystatin C Control Kit REF A52765.

Measuring patient samples

When a valid calibration has been performed and the control values are within the valid range, serum or plasma samples may be measured. Ensure that minimum volume of sample is present and assay the samples according to the instructions given in the Beckman Coulter® AU Systems instrument manuals.

Results

The results are calculated automatically by the Beckman Coulter® AU Systems. The results are presented in mg/L.

Relevant calculations

GFR prediction calculation

Several cystatin C based prediction equations for calculation of GFR for adults and children have been published. It should be noted that these formulas were evaluated with different cystatin C assays (Particle Enhanced Nephelometric Immunoassay PENIA or Particle Enhanced Turbidimetric Immunoassay PETIA) and may reveal inaccurate GFR results if an inappropriate combination of formula and assay is used. For calculation of GFR from cystatin C values measured with the Gentian assay the following prediction equation is recommended using mg/L as the unit factor [6]. The equation is valid for people above 14 years.

$$\text{GFR [mL/min/1.73 m}^2\text{]} = \frac{79.901}{\text{Cystatin C (mg/L)}^{1.4389}}$$

Limitations of the procedure

The materials should not be used past the expiration date.

Expected values

Gentian follows the CLSI Guideline, C28 [7] to determine the transferability of the reference interval. The reference interval is based on a reference interval study performed at Växjö Hospital, Sweden, including serum samples from 136 self-declared healthy subjects 20-80 years of age. The samples were analyzed for cystatin C on the AU2700 platform. The reference interval was calculated non-parametrically and was determined to be 0.53-1.01 mg/L. This represents the central 95 % of the whole population tested. It is recommended that every laboratory should determine a local reference interval since values may vary depending on the population tested.

Specific performance characteristics AU5800

All results refer to validation of the Gentian Cystatin C Immunoassay on an AU5800 instrument at one site with one lot of reagents, unless otherwise stated.

Measuring range

The measuring range of the Gentian Cystatin C Immunoassay was found to be 0.49-7.07 mg/L. The exact measuring range is specific to the calibrator, please refer to the analytical value sheet for the lot specific calibrator values available on www.gentian.com.

Interference

The interference study was designed in accordance with the protocol EP07 from CLSI [8]. Previously, no significant interference was detected with the drugs tested as recommended in a publication by Sonntag and Scholer [9]. There is no RF interference present in the Gentian Cystatin C Immunoassay because the antibodies are made using avian antibodies (chicken) [10].

Potential interferents	Concentration with no interference
Hemoglobin	6 g/L
Intralipid	10 g/L
Bilirubin	0.4 g/L

Precision

The Gentian Cystatin C Immunoassay was used in a study designed in accordance with CLSI protocol EP05 [11]. Three serum pools and 2 control levels were measured on the Beckman Coulter® AU5800 system (n=20).

Sample ID	Mean [mg/L]	Within run CV [%]	Total CV [%]
P1	0.90	0.82	1.96
P2	5.29	0.49	2.10
P3	2.08	0.43	1.62
CL	0.86	1.10	3.42
CH	2.91	0.81	2.40

Sensitivity

Using the Gentian Cystatin C Immunoassay on an AU5800 instrument, a lower limit of quantification was measured to be 0.23 mg/L.

Linearity

Using the Gentian Cystatin C Immunoassay, linearity was measured within acceptable limits in the range of 0.49-7.07 mg/L on the AU5800 system. Linearity samples above this range were not tested.

Hook effect

In a study on AU5800, the security zone for antigen excess extended up to 32 mg/L using the Gentian Cystatin C Immunoassay. No samples above this value were measured.

Analytical recovery

Using the Gentian Cystatin C Immunoassay on a Beckman Coulter® AU5800 instrument, a recovery of 96-100 % was observed.

Instrument variation

Instrument variation between Gentian Cystatin C Immunoassay on AU5800 and Architect c16000 instruments was measured and the results analyzed using Passing-Bablok regression analysis:

Passing-Bablok regression	N	Range of samples [mg/L]	Term	Coefficient
AU5800 vs. Architect	32	0.76 -1.88	Intercept	0.01
			Slope	0.95

Cystatin C

Specific performance characteristics AU680

All results refer to validation of the Gentian Cystatin C Immunoassay on an AU680 instrument at one site with one lot of reagents, unless otherwise stated.

Measuring range

The measuring range of the Gentian Cystatin C Immunoassay was found to be 0.44-7.30 mg/L. The exact measuring range is specific to the calibrator, please refer to the analytical value sheet for the lot specific calibrator values available on www.gentian.com.

Interference

The interference study was designed in accordance with the protocol EP07 from CLSI [8]. Previously, no significant interference was detected with the drugs tested as recommended in a publication by Sonntag and Scholer [9]. There is no RF interference present in the Gentian Cystatin C Immunoassay because the antibodies are made using avian antibodies (chicken) [10].

Potential interferents	Concentration with no interference
Hemoglobin	8.5 g/L
Intralipid	16 g/L
Bilirubin	0.2 g/L

Precision

The Gentian Cystatin C Immunoassay was used in study designed in accordance with CLSI protocol EP05 [11]. Four serum pools were measured on the Beckman Coulter® AU680 system (n=20).

Sample ID	Mean [mg/L]	Within run CV [%]	Between run CV [%]	Total CV [%]
P1	0.75	0.79	2.08	2.44
P2	1.96	0.43	1.73	1.88
P3	0.80	1.09	1.35	2.00
P4	4.98	0.67	1.00	1.57
CL	1.07	0.42	1.66	2.26
CH	3.28	0.25	1.00	1.51

Sensitivity

Using the Gentian Cystatin C Immunoassay on an AU680 instrument, a lower limit of quantification was measured to be 0.28 mg/L.

Linearity

Using the Gentian Cystatin C Immunoassay, linearity was measured within acceptable limits in the range of 0.44-7.30 mg/L on the AU680 system.

Hook effect

In a study on AU680, the security zone for antigen excess extended up to 12 mg/L using the Gentian Cystatin C Immunoassay.

Analytical recovery

Using the Gentian Cystatin C Immunoassay on a Beckman Coulter® AU680 instrument, a recovery of 86-92 % was observed.

Instrument variation

Instrument variation between Gentian Cystatin C Immunoassay on AU680 and Architect c16000 instruments was measured and the results analyzed using Passing-Bablok regression analysis:

Passing-Bablok regression	N	Range of samples [mg/L]	Term	Coefficient
AU680 Vs. Architect	40	0.71 – 6.38	Intercept	0.03
			Slope	0.95

Specific performance characteristics AU480

All results refer to validation of the Gentian Cystatin C Immunoassay on an AU480 instrument at one site with one lot of reagents, unless otherwise stated.

Measuring Range

The measuring range of the Gentian Cystatin C Immunoassay was found to be 0.43-7.32 mg/L. The exact measuring range is specific to the calibrator, please refer to the analytical value sheet for the lot specific calibrator values available on www.gentian.com.

Interference

The interference study was designed in accordance with the protocol EP07 from CLSI [8]. Previously, no significant interference was detected with the drugs tested as recommended in a publication by Sonntag and Scholer [9]. There is no RF interference present in the Gentian Cystatin C Immunoassay because the antibodies are made using avian antibodies (chicken) [10].

Potential interferents	Concentration with no interference
Hemoglobin	10 g/L
Intralipid	15 g/L
Bilirubin	0.6 g/L

Precision

The Gentian Cystatin C Immunoassay was used in a n study designed in accordance with CLSI protocol EP05 [11]. Three serum pools and 2 control levels were measured on the Beckman Coulter® AU480 system (n=12).

Sample ID	Mean [mg/L]	Within run CV [%]	Between run CV [%]	Total CV [%]
P1	1.09	1.57	1.21	3.60
P2	3.65	0.67	0.62	1.82
P3	1.24	1.73	0.00	3.47
CL	0.87	3.10	0.00	3.72
CH	3.39	1.18	0.94	3.03

Cystatin C

Sensitivity

Using the Gentian Cystatin C Immunoassay on an AU480 instrument, a lower limit of quantification was measured as 0.43 mg/L.

Linearity

Using the Gentian Cystatin C Immunoassay, linearity was measured within acceptable limits in the range of 0.40–7.32 mg/L on the AU480 system.

Hook effect

In a study on AU480, the security zone for antigen excess extended up to 9.4 mg/L using the Gentian Cystatin C Immunoassay.

Analytical recovery

Using the Gentian Cystatin C Immunoassay on a Beckman Coulter® AU480 instrument, a recovery of 90–95 % was observed.

Instrument variation

Instrument variation between Gentian Cystatin C Immunoassay on AU480 and Architect c16000 instruments was measured and the results analyzed using Passing-Bablok regression analysis:

Passing-Bablok regression	N	Range of samples (mg/L)	Term	Coefficient
AU480 vs. Architect	40	0.71 – 6.38	Intercept	0.03
			Slope	0.95

Performance characteristics DxC 500 AU/DxC 500i

All results refer to validation of the Gentian Cystatin C Immunoassay on a DxC 500 AU instrument at one site with one lot of reagents, unless otherwise stated.

Measuring range

The measuring range of the Gentian Cystatin C Immunoassay was found to be 0.38–7.84 mg/L. The exact measuring range is specific to the calibrator, please refer to the analytical value sheet for the lot specific calibrator values available on www.gentian.com.

Interference

The interference study was designed in accordance with the protocol EP07 from CLSI [8]. Previously, no significant interference was detected with the drugs tested as recommended in a publication by Sonntag and Scholer [9]. There is no RF interference present in the Gentian Cystatin C Immunoassay because the antibodies are made using avian antibodies (chicken) [10].

Potential interferents	Concentration with no interference
Haemoglobin	8 g/L
Intralipid	10 g/L
Bilirubin	0.2 g/L

Precision

The Gentian Cystatin C Immunoassay was used in a n study designed in accordance with CLSI protocol EP05 [11]. 3 serum pools and 2 controls were measured 2 times with 2 replicates (n=80).

Sample ID	Mean [mg/L]	Within run CV [%]	Between run CV [%]	Total CV [%]
P1	0.87	0.56	1.46	2.41
P2	1.60	0.80	1.63	2.43
P3	6.37	0.73	1.63	3.66
CL	1.00	0.68	0.61	2.00
CH	3.48	0.46	0.55	1.57

Sensitivity

Using the Gentian Cystatin C Immunoassay on a DxC 500 AU instrument, a lower limit of quantification was measured as 0.32 mg/L. The study was designed in accordance with EP17 [12].

Linearity

Using the Gentian Cystatin C Immunoassay, linearity was measured within acceptable limits in the range of 0.38–7.84 mg/L on the DxC 500 AU system. Linearity samples above this range were not tested.

Hook effect

In a study on DxC 500 AU, the security zone for antigen excess extended up to 25.7 mg/L using the Gentian Cystatin C Immunoassay.

Analytical recovery

Using the Gentian Cystatin C Immunoassay on a DxC 500 AU instrument, a recovery of 102–109 % was observed.

Instrument variation

Instrument variation between Gentian Cystatin C Immunoassay on the DxC 500 AU and AU5800 instruments was measured and the results analyzed using Passing-Bablok regression analysis:

n	Range of samples [mg/L]	Term	Coefficient	95% CI
42	0.57– 5.72	Intercept	-0.01	[-0.05, 0.03]
		Slope	1.00	[0.97, 1.04]

Performance characteristics DxC 700 AU

All results refer to validation of the Gentian Cystatin C Immunoassay on a DxC 700 AU instrument at one site with one lot of reagents, unless otherwise stated.

Measuring range

The measuring range of the Gentian Cystatin C Immunoassay was found to be 0.40–8.07 mg/L. The exact measuring range is specific to the calibrator, please refer to the analytical value sheet for the lot specific calibrator values available on www.gentian.com.

Interference

The interference study was designed in accordance with the protocol EP07 from CLSI [8]. Previously, no significant interference was detected with the drugs tested as recommended in a publication by Sonntag and Scholer [9]. There is no RF interference present in the Gentian Cystatin C Immunoassay because the antibodies are made using avian antibodies (chicken) [10].

Potential interferents	Concentration with no interference
Hemoglobin	10 g/L
Intralipid	10 g/L
Bilirubin	0.2 g/L

Cystatin C



Precision

The Gentian Cystatin C Immunoassay was used in a study designed in accordance with CLSI protocol EP05 [11]. Three serum pools and 2 control levels were measured on the Beckman Coulter® DxC 700 AU system (n=80).

Sample ID	Mean [mg/L]	Within run CV [%]	Between Run CV [%]	Total CV [%]
P1	0.73	0.58	0.00	0.75
P2	1.70	0.49	0.28	0.59
P3	6.13	0.44	0.18	0.60
CL	0.91	0.67	0.60	1.04
CH	3.44	0.39	0.81	0.90

Sensitivity

Using the Gentian Cystatin C Immunoassay on an DxC 700 AU instrument, a lower limit of quantification was measured as 0.40 mg/L. The study was designed in accordance with EP17 [12].

Linearity

Using the Gentian Cystatin C Immunoassay, linearity was measured within acceptable limits in the range of 0.40-8.07 mg/L on the DxC 700 AU system. Linearity samples above this range were not tested.

Hook effect

In a study on DxC 700 AU, the security zone for antigen excess extended up to 10 mg/L using the Gentian Cystatin C Immunoassay.

Analytical recovery

Using the Gentian Cystatin C Immunoassay on a Beckman Coulter® DxC 700 AU instrument, a recovery of 104 – 105 % was observed.

Instrument variation

Instrument variation between Gentian Cystatin C Immunoassay on DxC 700 AU and Architect c4000, and between DxC 700 AU and AU5800 instruments was measured and the results analyzed using Passing-Bablok regression analysis:

Instrument	N	Range of samples [mg/L]	Term	Co-efficient
Architect	40	0.60-6.27	Intercept	0.02
			Slope	0.96
AU 5800	40	0.59-6.22	Intercept	0.00
			Slope	1.00

Bibliography

1. Abrahamson M et al: Biochem J 1990;268:287-94
2. Laterza OF et al: Clin Chem 2002;48:63-99.
3. Grubb AO. Adv Clin Chem 2000;35:63-99.
4. Filler G et al: Clin Biochem 2005 ;38 :1-8
5. Shlipak M.G, et al: Clinical Chemistry 57: 737-745, 2011
6. Flodin M et al: Scand J Clin Lab Invest 2007;67:560-567
7. CLSI. *Defining, Establishing, and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline – Third Edition*. CLSI document C28-A3c. Wayne, PA: Clinical Laboratory Standards Institute; 2008.
8. CLSI. *Interference Testing in Clinical Chemistry*. 3rd ed. CLSI guideline EP07. Wayne, PA: Clinical Laboratory Standards Institute; 2018
9. Sonntag O, Scholer A. Ann Clin Biochem 2001;38:376-85.
10. Larsson A et al: J Immunol Methods. 1988 Apr 6;108(1-2):205-8
11. CLSI. *Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline – Third Edition*. CLSI document EP05-A3. Wayne, PA: Clinical Laboratory Standards Institute; 2014.

12. CLSI. *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline – Second Edition*. CLSI document EP17-A2. Wayne, PA: Clinical and Laboratory Standards Institute ;2012.



Gentian AS

Bjornasveien 5
N-1596 Moss
Norway
TEL: +47 99 33 99 05
www.gentian.com



RxOnly

Date of issue

21 MAR. 2025

Additional information

For more detailed information on AU Systems, refer to the appropriate system manual. Since Beckman Coulter® does not manufacture the reagent or perform quality control or other tests on individual lots, Beckman Coulter® cannot be held responsible for the quality of the data obtained which is caused by performance of the reagent, any variation between lots of reagent, or protocol changes by the manufacturer.


Serious incidents

Please notify the distributor and your competent authority if any serious incidents have occurred in relation to the device.

Shipping damage

Please notify your distributor if this product is received damaged. For technical assistance please contact your local Beckman Coulter® representative.

Symbols key



Temperature limit



Use by date



Consult instructions for use



Manufacturer



CE mark with Notified Body number



UKCA mark



In Vitro Diagnostic medical device



Lot number



Catalogue number



Unique Device Identifier



Contents

Cystatin C



R1

R1 Assay Buffer

R2

R2 Immunoparticles

RxOnly

Caution: Federal law restricts this device to sale by or on the order of a physician.

Cystatin C



Instrument settings for the Gentian Cystatin C Immunoassay

Cystatin C AU5800 application settings

System Reagent: B08179 Reagent ID: 228

Parameters		Specific Test Parameters	
General	LIH	ISE	Calculated Test Range
Test Name: CysC < > Type: Serum Operation: Yes			
Sample Volume	2 μL	Dilution	0 μL OD Limit
Pre-Dilution Rate	1	Diluent Bottle	Outside
Rgt. Volume	150 μL	Dilution	0 μL Reagent OD Limit
R1-2	150 μL	Dilution	0 μL First Low: -2.0 High: 2.0
R2(R2-1)	30 μL	Dilution	10 μL Last Low: 0.49 High: 7.07
Common Rgt. Type		Name	
Wavelength	540 nm	Sec.	
Method	End Point	Correlation Factor A	1.00
Reaction Slope	13	Onboard Stability Period	60** Day
Measuring Point1 First	13	LIH Influence Check	
Measuring Point2 First	27	Lipemia	
Linearity Limit		Icterus	
Lag Time Check		Hemolysis	

Parameters		Specific Test Parameters	
General	LIH	ISE	Calculated Test Range
Test Name: CysC < > Type: Serum			
Value/Flag: #			
Level: Low # High #			
Specific Ranges: From To			
1.	Sex	Year	Month
2.	Sex	Year	Month
3.	Sex	Year	Month
4.	Sex	Year	Month
5.	Sex	Year	Month
6.	Sex	Year	Month
7.	Standard demographics		
8.	Not within expected values		
Panic Value	Low	High	Unit: mg/L Decimal Places: #

Parameters		Calibration Parameters	
Calibrators	Calibration Specific	STAT Table Calibration	
General	ISE		
Test Name: CysC < > Type: Serum Cuvette:			
Calibration Type: 6AB Formula: Spline Counts: #			
<Calibrator Parameters>			
Point 1:	1	OD	*
Point 2:	2	OD	*
Point 3:	3	OD	*
Point 4:	4	OD	*
Point 5:	5	OD	*
Point 6:	6	OD	*
Point 7:	7	OD	*
Point 8:	8	OD	*
Point 9:	9	OD	*
Point 10:	10	OD	*
Slope Check: +			
Allowance Range Check			
☐ Reagent Blank			
☐ Calibration			
Advanced Calibration Operation: #			
Interval (RB/ACAL): #			
<Point Cal. For Master Curve>			
No. of Correction Points: Use Master Curve:			
☐ Lot Calibration			
Point-1	Calibrator	OD	Conc
Point-2	Calibrator	OD	Conc
Stability: Reagent Blank: 28*** Day: 0 Hour			
Calibration: 28*** Day: 0 Hour			
MB Type Factor: 1-Point Calibration Point with Conc-0			

User defined

* Lot specific, see analytical value sheet available on www.gentian.com

** based on results from instrument AU400 (Beckman Coulter®)

*** based on results from instrument DxC 700 AU

Cystatin C



Cystatin C AU680 application settings

System Reagent: B08179

Reagent ID: 228

Parameters		Specific Test Parameters	
General	LIH	ISE	Calculated Test Range
Test Name: CysC < > Type: Serum Operation: Yes			
Sample Volume: 2 μL Dilution: 0 μL OD Limit: Pre-Dilution Rate: 1 μL Dilution: 0 μL Reagent OD Limit: Rgt. Volume: 150 μL Dilution: 0 μL First: -2.0 Low: 2.0 High: R2(R2-1): 30 μL Dilution: 10 μL Last: 0.44 Low: 7.30 High: Common Rgt. Type: Name: Dynamic Range Low: 0.44 High: 7.30 Wavelength: Pri: 540 nm Sec: nm Correlation Factor A: 1.00 B: 0.00 Method: End Point Factor for Maker A: Reaction Slope: + Onboard Stability Period: 60** Day: Hour: Measuring Point1 First: 13 Last: 27 LIH Influence Check: Measuring Point2 First: Last: Lipemia: Linearity Limit: % Icterus: Lag Time Check: Hemolysis:			

Parameters		Specific Test Parameters																																																																									
General	LIH	ISE	Calculated Test Range																																																																								
Test Name: CysC < > Type: Serum																																																																											
Value/Flag: # Low: High: Level: Specific Ranges: From To <table><thead><tr><th></th><th>Sex</th><th>Year</th><th>Month</th><th>Year</th><th>Month</th><th>Low</th><th>High</th></tr></thead><tbody><tr><td>☐ 1.</td><td> </td><td> </td><td> </td><td> </td><td> </td><td> </td><td> </td></tr><tr><td>☐ 2.</td><td> </td><td> </td><td> </td><td> </td><td> </td><td> </td><td> </td></tr><tr><td>☐ 3.</td><td> </td><td> </td><td> </td><td> </td><td> </td><td> </td><td> </td></tr><tr><td>☐ 4.</td><td> </td><td> </td><td> </td><td> </td><td> </td><td> </td><td> </td></tr><tr><td>☐ 5.</td><td> </td><td> </td><td> </td><td> </td><td> </td><td> </td><td> </td></tr><tr><td>☐ 6.</td><td> </td><td> </td><td> </td><td> </td><td> </td><td> </td><td> </td></tr><tr><td>☐ 7.</td><td colspan="6">No demographics</td><td> </td></tr><tr><td>☐ 8.</td><td colspan="6">Not within expected values</td><td> </td></tr></tbody></table> Panic Value: Low: High: Unit: mg/L Decimal Places: #					Sex	Year	Month	Year	Month	Low	High	☐ 1.								☐ 2.								☐ 3.								☐ 4.								☐ 5.								☐ 6.								☐ 7.	No demographics							☐ 8.	Not within expected values						
	Sex	Year	Month	Year	Month	Low	High																																																																				
☐ 1.																																																																											
☐ 2.																																																																											
☐ 3.																																																																											
☐ 4.																																																																											
☐ 5.																																																																											
☐ 6.																																																																											
☐ 7.	No demographics																																																																										
☐ 8.	Not within expected values																																																																										

Parameters		Calibration Parameters																																																																																					
Calibrators	Calibration Specific	STAT Table Calibration																																																																																					
Test Name: CysC < > Type: Serum Use Serum Cal: ☐																																																																																							
Calibration Type: 6AB Formula: Spline Counts: # <Calibrator Parameters> <table><thead><tr><th>Calibrator</th><th>OD</th><th>Conc</th><th>Factor Range</th><th>Low</th><th>High</th></tr></thead><tbody><tr><td>Point 1:</td><td> </td><td> </td><td> </td><td> </td><td> </td></tr><tr><td>Point 2:</td><td> </td><td> </td><td> </td><td> </td><td> </td></tr><tr><td>Point 3:</td><td> </td><td> </td><td> </td><td> </td><td> </td></tr><tr><td>Point 4:</td><td> </td><td> </td><td> </td><td> </td><td> </td></tr><tr><td>Point 5:</td><td> </td><td> </td><td> </td><td> </td><td> </td></tr><tr><td>Point 6:</td><td> </td><td> </td><td> </td><td> </td><td> </td></tr><tr><td>Point 7:</td><td> </td><td> </td><td> </td><td> </td><td> </td></tr><tr><td>Point 8:</td><td> </td><td> </td><td> </td><td> </td><td> </td></tr><tr><td>Point 9:</td><td> </td><td> </td><td> </td><td> </td><td> </td></tr><tr><td>Point 10:</td><td> </td><td> </td><td> </td><td> </td><td> </td></tr></tbody></table> Slope Check: + Allowance Range Check: ☐ Reagent Blank: Calibration: Advanced Calibration Operation: Interval (RB/ACAL): <Point Cal. For Master Curve> No. of Correction Points: Use Master Curve: ☐ Lot Calibration: ☐ <table><thead><tr><th>Calibrator</th><th>OD</th><th>Conc</th><th>OD Range</th><th>Low</th><th>High</th></tr></thead><tbody><tr><td>Point-1</td><td> </td><td> </td><td> </td><td> </td><td> </td></tr><tr><td>Point-2</td><td> </td><td> </td><td> </td><td> </td><td> </td></tr></tbody></table> Stability: Reagent Blank: 28*** Day: Hour: Calibration: 28*** Day: Hour: MB Type Factor: 1-Point Calibration Point: ☐ with Conc-0: ☐				Calibrator	OD	Conc	Factor Range	Low	High	Point 1:						Point 2:						Point 3:						Point 4:						Point 5:						Point 6:						Point 7:						Point 8:						Point 9:						Point 10:						Calibrator	OD	Conc	OD Range	Low	High	Point-1						Point-2					
Calibrator	OD	Conc	Factor Range	Low	High																																																																																		
Point 1:																																																																																							
Point 2:																																																																																							
Point 3:																																																																																							
Point 4:																																																																																							
Point 5:																																																																																							
Point 6:																																																																																							
Point 7:																																																																																							
Point 8:																																																																																							
Point 9:																																																																																							
Point 10:																																																																																							
Calibrator	OD	Conc	OD Range	Low	High																																																																																		
Point-1																																																																																							
Point-2																																																																																							

User defined

* Lot specific, see analytical value sheet available on www.gentian.com

** based on results from instrument AU400 (Beckman Coulter®)

*** based on results from instrument Dx C 700 AU

Cystatin C



Cystatin C AU480 application settings

System Reagent: B08179

Reagent ID: 228

Parameters		Specific Test Parameters																																																																																																								
General	LIH	ISE	Calculated Test	Range																																																																																																						
Test Name: CysC < > Type: Serum Operation: Yes																																																																																																										
<table><tr><td>Sample Volume</td><td> 2 μL</td><td>Dilution</td><td> 0 μL</td><td>OD Limit</td><td></td></tr><tr><td>Pre-Dilution Rate</td><td> 1 </td><td></td><td></td><td>Min.OD</td><td> </td></tr><tr><td>Rgt. Volume</td><td> 150 μL</td><td>Dilution</td><td> 0 μL</td><td>Reagent OD Limit</td><td> </td></tr><tr><td></td><td></td><td></td><td></td><td>First</td><td> -2.0 </td></tr><tr><td></td><td></td><td></td><td></td><td>Last</td><td> 2.0 </td></tr><tr><td></td><td></td><td></td><td></td><td>Low</td><td> </td></tr><tr><td></td><td></td><td></td><td></td><td>High</td><td> </td></tr><tr><td>R2(R2-1)</td><td> 30 μL</td><td>Dilution</td><td> 10 μL</td><td>Dynamic Range Low</td><td> 0.43 </td></tr><tr><td></td><td></td><td></td><td></td><td>High</td><td> 7.32 </td></tr><tr><td>Common Rgt. Type</td><td> </td><td>Name</td><td> </td><td>Correlation Factor A</td><td> 1.00 </td></tr><tr><td>Wavelength</td><td> 540 nm</td><td>Sec.</td><td> </td><td>Factor for Maker A</td><td> 0.00 </td></tr><tr><td>Method</td><td> End Point </td><td></td><td></td><td></td><td></td></tr><tr><td>Reaction Slope</td><td> + </td><td></td><td></td><td>Onboard Stability Period</td><td> 60** Day</td></tr><tr><td>Measuring Point1 First</td><td> 13 </td><td>Last</td><td> 27 </td><td>LIH Influence Check</td><td> </td></tr><tr><td>Measuring Point2 First</td><td> </td><td>Last</td><td> </td><td>Lipemia</td><td> </td></tr><tr><td>Linearity Limit</td><td> </td><td></td><td></td><td>Icterus</td><td> </td></tr><tr><td>Lag Time Check</td><td> </td><td></td><td></td><td>Hemolysis</td><td> </td></tr></table>					Sample Volume	2 μL	Dilution	0 μL	OD Limit		Pre-Dilution Rate	1			Min.OD		Rgt. Volume	150 μL	Dilution	0 μL	Reagent OD Limit						First	-2.0					Last	2.0					Low						High		R2(R2-1)	30 μL	Dilution	10 μL	Dynamic Range Low	0.43					High	7.32	Common Rgt. Type		Name		Correlation Factor A	1.00	Wavelength	540 nm	Sec.		Factor for Maker A	0.00	Method	End Point					Reaction Slope	+			Onboard Stability Period	60** Day	Measuring Point1 First	13	Last	27	LIH Influence Check		Measuring Point2 First		Last		Lipemia		Linearity Limit				Icterus		Lag Time Check				Hemolysis	
Sample Volume	2 μL	Dilution	0 μL	OD Limit																																																																																																						
Pre-Dilution Rate	1			Min.OD																																																																																																						
Rgt. Volume	150 μL	Dilution	0 μL	Reagent OD Limit																																																																																																						
				First	-2.0																																																																																																					
				Last	2.0																																																																																																					
				Low																																																																																																						
				High																																																																																																						
R2(R2-1)	30 μL	Dilution	10 μL	Dynamic Range Low	0.43																																																																																																					
				High	7.32																																																																																																					
Common Rgt. Type		Name		Correlation Factor A	1.00																																																																																																					
Wavelength	540 nm	Sec.		Factor for Maker A	0.00																																																																																																					
Method	End Point																																																																																																									
Reaction Slope	+			Onboard Stability Period	60** Day																																																																																																					
Measuring Point1 First	13	Last	27	LIH Influence Check																																																																																																						
Measuring Point2 First		Last		Lipemia																																																																																																						
Linearity Limit				Icterus																																																																																																						
Lag Time Check				Hemolysis																																																																																																						
Parameters																																																																																																										
Specific Test Parameters																																																																																																										
General																																																																																																										
LIH																																																																																																										
ISE																																																																																																										
Calculated Test																																																																																																										
Range																																																																																																										
Test Name: CysC < > Type: Serum																																																																																																										
Value/Flag: # Low # High #																																																																																																										
Level																																																																																																										
Specific Ranges: From # To #																																																																																																										
Panic Value																																																																																																										
Low # High #																																																																																																										
1. # # # # # # # #																																																																																																										
2. # # # # # # # #																																																																																																										
3. # # # # # # # #																																																																																																										
4. # # # # # # # #																																																																																																										
5. # # # # # # # #																																																																																																										
6. # # # # # # # #																																																																																																										
7. No demographics																																																																																																										
8. Not within expected values																																																																																																										
Unit mg/L Decimal Places #																																																																																																										
Parameters																																																																																																										
Calibration Parameters																																																																																																										
Calibrators																																																																																																										
Calibration Specific																																																																																																										
STAT Table Calibration																																																																																																										
General																																																																																																										
ISE																																																																																																										
Test Name: CysC < > Type: Serum Use Serum Cal. ☐																																																																																																										
Calibration Type: 6AB Formula: Spline Counts: #																																																																																																										
<Calibrator Parameters>																																																																																																										
Calibrator																																																																																																										
OD																																																																																																										
Conc																																																																																																										
Factor Range																																																																																																										
Low																																																																																																										
High																																																																																																										
Slope Check																																																																																																										
Allowance Range Check																																																																																																										
Point 1: 1 OD Conc Factor Range																																																																																																										
Point 2: 2 OD Conc Factor Range																																																																																																										
Point 3: 3 OD Conc Factor Range																																																																																																										
Point 4: 4 OD Conc Factor Range																																																																																																										
Point 5: 5 OD Conc Factor Range																																																																																																										
Point 6: 6 OD Conc Factor Range																																																																																																										
Point 7: 7 OD Conc Factor Range																																																																																																										
Point 8: 8 OD Conc Factor Range																																																																																																										
Point 9: 9 OD Conc Factor Range																																																																																																										
Point 10: 10 OD Conc Factor Range																																																																																																										
<Point Cal. For																																																																																																										
No. of Correction Points																																																																																																										
Use Master Curve																																																																																																										
Lot Calibration																																																																																																										
Master Curve																																																																																																										
Calibrator																																																																																																										
OD																																																																																																										
Conc																																																																																																										
Low																																																																																																										
High																																																																																																										
Stability																																																																																																										
Reagent Blank																																																																																																										
Calibration																																																																																																										
Interval (RB/ACAL)																																																																																																										
MB Type Factor: 1-Point Calibration Point																																																																																																										

User defined

* Lot specific, see analytical value sheet available on www.gentian.com

** based on results from instrument AU400 (Beckman Coulter®)

*** based on results from instrument DxH 700 AU

Cystatin C



Cystatin C DxC 500 AU/DxC 500i application settings

System Reagent: B08179

Reagent ID: 228

TEST CONFIGURATION & CHEMISTRY DETAILS									
Assay Name	Test	Rev	Discipline				Chemistry		
Test ID	CYS		Calculated Result				☐		
LIS Code	CYS								
UNITS AND RANGE SETTINGS									
Use Settings from	Serum	Units	mg/L	Decimal Places	X.XX	Plasma			
Test Kind	General	Revision	01	☒ Multi Reagent Switch					
Reagent Name	CYS	Reagent ID	228	☐ FSE Test					
ABB Name	CYS1G	Parameter Long Name	Cystatin C B08179 CYS1G CYSC Serum						
Region	☒ US ☒ OUS ☒ AP ☐ JP ☒ EU ☐ Other								
GENERAL PARAMETERS									
SAMPLE VOLUME			REACTION OD LIMIT						
Sample Volume	2.0	µL	Dilution	0	µL	Low		High	
Predilution Rate	1								
REAGENT VOLUME			REACTION BLANK OD LIMIT						
R1-1	150	µL	Dilution	0	µL	First: Low	-2.0000	High	2.0000
R2-1	30	µL	Dilution	10	µL	Last: Low	-2.0000	High	2.0000
WAVELENGTH			ANALYTICAL MEASURING RANGE						
Primary	540	nm	Secondary	None	nm	Low	0.38	High	7.84
METHOD			MANUFACTURER FACTOR						
END			A			1	B	0	
REACTION SLOPE			REAGENT ONBOARD STABILITY						
+			60**			Days	0**	Hours	
MEASURING POINT			LIH INFLUENCE CHECK						
Point 1: First	13	Last	27	☐ Perform LIH check					
Point 2: First		Last		Lipemia	+				
Linearity Limit			Icterus			+			
			Hemolysis			+			
Lag Time Check			☐ Perform Lag Time Check						
CALIBRATION PARAMETERS									
Base Unit	Decimal Place	Unit 1	Factor 1	Unit 2	Factor 2	Unit 3	Factor 3	Unit 4	Factor 4
mg/L	2	None	0	None	0	None	0	None	0
CALIBRATOR SPECIFIC									
Calibration Type			6AB	Counts	2	☐ Use highest calibrator for Upper AMR			
Formula			Spline	MB Factor					
Calibrator Name				Positive Cutoff					
Add			CYS						
☒ SLOPE CHECK			Number of Levels	6					
Slope Check			+						
STABILITY AND INTERVAL									
Reagent Blank Stability	28	Days	0	Hours	Interval	Lot			
Calibration Stability	28	Days	0	Hours	Interval	Lot			
OD DELTA CHECK									
☐ Reagent Blank			0.0000						
☐ Calibration			0.0000						
PROZONE CHECK PARAMETERS									
☐ Logic Check 1			☐ Logic Check 2			☐ Logic Check 3			
Check Points	Decision Values	Check Points	Decision Values	Check Points	Decision Values	Check Points	Decision Values	Check Points	Decision Values
Point 1	#	Point 1	#	Point 1	#	Point 1	#	Point 1	#
Point 2	#	Point 2	#	Point 2	#	Point 2	#	Point 2	#
Point 3	#	Point 3	#	Point 3	#	Point 3	#	Point 3	#
Limit Points	Limit 1	Limit Points	Limit 1	Limit Points	Limit 1	Limit Points	Limit 1	Limit Points	Limit 1
Limit 1	#	Limit 1	#	Limit 1	#	Limit 1	#	Limit 1	#
Limit 2	#	Limit 2	#	Limit 2	#	Limit 2	#	Limit 2	#
Check Pattern	Pattern	Check Pattern	Pattern	Check Pattern	Pattern	Check Pattern	Pattern	Check Pattern	Pattern
Pattern	#	Pattern	#	Pattern	#	Pattern	#	Pattern	#

User defined

* Lot specific, see analytical value sheet available on www.gentian.com

** Based on results from instrument AU400 (Beckman Coulter®)

Cystatin C



Cystatin C Dx C 700 AU application settings

System Reagent: B08179 Reagent ID: 228

General		LIH		ISE		Calculated Test		Range	
Test Name: CYS1G		Test No		Type: Serum		Operation: Yes			
Sample Volume: 2.0 µL		Dilution: 0 µL		OD Limit		Min. OD		Max OD	
Pre-Dilution Rate: 1						Reagent OD Limit 1st		Low: -2.0000 High: 2.0000	
Reagent Volume R1 (R1-1): 150 µL		Dilution: 0 µL				Last		Low: -2.0000 High: 2.0000	
R1-2: 30 µL		Dilution: 10 µL				Analytical Measuring Range		Low: 0.40 High: 8.07	
R2 (R2-1): 30 µL						Correlation Factor		A: 1 B: 0	
Common Reagent Type: None		Name: None				Manufacturer Factor		A: 1 B: 0	
Wavelength Pri: 540 nm		Sec: None nm				Onboard Stability Period: 60** Day		0 Hour	
Method: END						LIH Influence Check: No			
Reaction Slope: +						Lipemia			
Measuring Point-1 1st: 13		Last: 27				Icterus			
Measuring Point-2 1st: 13		Last: 27				Hemolysis			
Linearity Limit: %									
Lag Time Check: %									

General		LIH		ISE		Calculated Test		Range	
Test Name: CYS1G		Test No		Type: Serum					
Value/Flag: Value		Level		Low: -99999.99		High: 99999.99			
Specific Ranges									
Sex: #		Year: #		Month: #		Other Type: None		Low: # High: #	
□ 1: #		#		#		None		# #	
□ 2: #		#		#		None		# #	
□ 3: #		#		#		None		# #	
□ 4: #		#		#		None		# #	
□ 5: #		#		#		None		# #	
□ 6: #		#		#		None		# #	
7: Standard demographics								# #	
8: Not within expected values								# #	
Critical Limits		Low: #		High: #		Unit: mg/L		Select Decimal Places: 2	

Calibrators		General		ISE	
Test Name: CYS1G				Type: Serum	
		□ Use Serum Cal.			
Calibration Type: 6AB		Formula: Spline		Counts: 2	
<Calibrator Parameters>				Slope Check: +	
				Allowable Range Check	
				□ Reagent Blank	
				□ Calibration	
				Advanced Calibration	
				Operation: No	
				Interval (RB):	
				Interval (ACAL):	
				Stability	
				Reagent Blank: 28 Day 0 Hour	
				Calibration: 28 Day 0 Hour	

Point	Calibrator	OD	Conc	Range	
				Low	High
Point-1	CYSC Calibrator Level 1		*	-2.0000	2.0000
Point-2	CYSC Calibrator Level 2		*	-2.0000	2.0000
Point-3	CYSC Calibrator Level 3		*	-2.0000	2.0000
Point-4	CYSC Calibrator Level 4		*	-2.0000	2.0000
Point-5	CYSC Calibrator Level 5		*	-2.0000	2.0000
Point-6	CYSC Calibrator Level 6		*	-2.0000	2.0000
Point-7					

MB Type Factor: 1-Point Calibration Point: None □ with Conc-0

User defined

* Lot specific, see analytical value sheet available on www.gentian.com

** based on results from instrument AU400 (Beckman Coulter®)